



SUNMARK
INSTRUMENTS

Mark Of Quality



Our Range of Best Selling Instruments

DENTAL INSTRUMENTS CATALOGUE

INSTRUCTIONS FOR REPROCESSING OF REUSABLE SURGICAL INSTRUMENTS

The following instructions are for all reusable medical Devices supplied by Sunmark Instruments Ltd, unless stated otherwise with the packaging of the product. These instructions are intended for use only by persons with the required specialist knowledge and training.

Warnings

- Follow instructions and warnings as issued by the manufacturers of any Decontaminants, disinfectants and cleaning agents. Wherever possible avoid use of mineral aci37 and harsh, abrasive agents.
- No part of the process shall exceed 140 °C.
- Some sensitive materials (e.g. Aluminium) are damaged by high alkaline solutions (pH>10).
- Devices with long, narrow cannula, hinges and blind holes require particular attention during cleaning.

Note: When reprocessing medical Devices, always handle with care, wearing protective clothing, gloves and eye-wear in accordance with local health & safety procedures.

Limitations on Reprocessing

- Repeated processing has minimal effect on these instruments.
- End of life normally Determined by wear and damage in use.
- Any specific limitations on the number of reprocessing cycles shall be ma33 available with instrument.

From Point of Use

Whenever possible, 36 not allow blood, 33bris or bodily flui37 to dry on instruments. For best results and to prolong the life of the medical 33vice reprocess immediately after use. If they cannot be reprocessed immediately, use an enzymatic foam spray cleaner to help prevent soil from drying.

Preparation for Decontamination

- Process all the instruments as soon as it is reasonably practical following use.
- Disassemble only where inten33d, without the use of tools unless specifically provi33d by the manufacturer. Instructions for dis assembly are available with the 33vice, where required.

Cleaning: Automated

- Use only either CE marked or validated washer-disinfector machines and low foaming, non-ionizing cleaning agents and Detergents following the manufacturer's instructions for use, warnings, concentrations and recommenDed cycles.
- Load instruments carefully, with any box joints and hinges open and so that any fenestration in instruments can drain.
- Place heavy instruments with care in the bottom of containers, taking care not to overload wash baskets.
- Place instruments with concave surfaces facing down to prevent pooling of water.
- Where possible, use appropriate attachments to flush insi33 reamers and 33vices with lumens or cannula.
- Ensure that soft, high purity water which is controlled for bacterial en36-toxins is used in the final rinse stage.

Note: Automated cleaning may not be suitable for all lumens and cannula in which case clean manually with a water jet gun, if available, and an appropriate brush (and stilette if proviDed) that reaches the Depth of the feature. After manually cleaning, pass all Devices through an automatic cleaning cycle to achieve disinfection.

Note: These instruments have been validated using a washer-disinfector cycle validated to incluDe two cold rinses at <35 °C, a Detergent cycle and a rinse cycle both at >50°C, a disinfection cycle operating at a temperature of between 80 °C and 90 °C for a minimum holding time of 1 minute (actual holding time in excess of 2 minutes 50 secon37) and a 20 minute drying cycle. The Detergent used was low foaming, non-ionizing spray wash Detergent cleaner (max 12pH) and the rinse aid a neutral pH low foaming, non-ionic surfactant with isopropyl alcohol.

Cleaning: Manual

- Manual cleaning is not a38ised if an automatic washer-disinfector is available. If the equipment is not available, use the following process:-
1. Use a double sink system (wash/rinse) Dedicated for instrument cleaning (not used for hand washing). Ensure that the water temperature does not exceed 35 of the re-processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the Desired result. This requires validation and routine monitoring of the process. Likewise any Deviation by the re-processor from the instructions proviDed must be properly evaluated for effectiveness and potential a38erse consequences.
 2. In the first sink, keeping the instrument submerged, with an autoclavable brush, apply CE marked cleaning solution to all surfaces until all soil has been removed. Pay particular attention to serrations, teeth, ratchets and hinges, always brushing away from the body and avoiding splashing. Ensure rongeurs and hinged instruments are thoroughly cleaned in both open and closed positions.
 3. In the second sink, rinse instruments thoroughly with soft, high purity water which is controlled for bacterial endotoxins, so that the water reaches all parts of the instrument, then carefully hand dry or use a drying cabinet.

Note: Manual Cleaning is NOT a disinfection process: When manual cleaning is used it may not be possible to disinfect the Device prior to further handling.

Cleaning Inspection

- After cleaning visually inspect, visually inspect all surfaces, cannulations, ratchets, joints, holes and lumens for complete removal of soil and flui37. If ANY soil or fluid is still visible, return the instrument for repeat Decontamination.

Maintenance

- Apply surgical grade lubricant to hinges, joints and moving parts as per the lubricant manufacturer's instructions.

Continued on next page

Inspection and Functional Testing

- Visually inspect and check: - all instruments for damage and wear; cutting edges are free from nicks and present a continuous edge; jaws and teeth align correctly; all articulated instruments have a smooth movement without excess play; locking mechanisms (such as ratchets) fasten securely and close easily; long, slender instruments are not distorted; any component parts fit and assemble correctly with mating components.
- Remove for repair and replacement any blunt, worn out, flaking, fractured or damaged instruments.

Note: If an instrument is returned to the manufacturer / supplier, the instrument must be Decontaminated and sterilized and be accompanied with the relevant documented evidence.

Packaging

- All instruments to be packed following local protocol in accordance with ISO 11607-1 or AAMI/CSR technique.

Sterilization (USA)

- Autoclave should comply with the requirements of, and be validated and maintained in accordance with ANSI/AAMI ST79.
- Pre-Vacuum moist heat sterilization operating at temperature 132 °C bar for a minimum holding time of 4 minutes always following the instructions of the machine manufacturer.
- When sterilizing multiple instruments in one autoclave cycle, ensure that the sterilizer manufacturers stated maximum load is not exceeded.
- Ensure instruments are dry before sterilization.

Sterilization (Outside USA)

- Autoclave should comply with the requirements of, and be validated and maintained in accordance with EN 285, EN 13060, EN ISO 17665.
- Either CE marked or validated vacuum autoclave operating at 134-137 °C bar for a minimum holding time of 3 minutes always following the instructions of the machine manufacturer.
- When sterilizing multiple instruments in one autoclave cycle, ensure that the sterilizer manufacturers stated maximum load is not exceeded.

Storage

- Ensure instruments are dry before storage, and stored in dry, clean conditions at an ambient room temperature.

Additional Information

- Other forms of cleaning (i.e. ultrasonic) and sterilization (i.e. Low temperature steam and Formaldehyde, Ethylene Oxide and Gas Plasma) are available. However, always follow the instructions for use as issued by the manufacturer and always consult with them if in any doubt over the suitability of any process used.
- Cleaning and sterilizing guidelines are available in HTM 2030 and HTM 2010. Contact: The NHS Estates Stationers Office Publications Cent for 33tails at www.tsonline.gov.uk. For further information contact: NHS Estates Information Center, Department of Health, 1 Trevelyan Square, Boar Lane, Lee37, LS1 6AE, UK or visit www.nhsstates.gov.uk.

Note: It is the responsibility of the re-processor to ensure that the reprocessing as actually performed using equipment, materials and personnel in the reprocessing facility achieve the Desired result. This requires validation and routine monitoring of the process. Likewise any Deviation by the re-processor from the instructions provided must be properly evaluated for effectiveness and potential adverse consequences.



MEASUREMENTS, STANDAR37, AND CONVERSIONS

Millimeter Equivalents

MM = INCHES	MM = INCHES	MM = INCHES
.1 = .00384	29.0 = 1.14173	65.0 = 2.55905
.2 = .00787	30.0 = 1.18110	66.0 = 2.59842
.3 = .01181	31.0 = 1.22047	67.0 = 2.63779
.4 = .01575	32.0 = 1.25984	68.0 = 2.67716
.5 = .01968	De.0 = 1.29921	69.0 = 2.71653
.6 = .02do2	34.0 = 1.De858	70.0 = 2.75590
.7 = .02756	35.0 = 1.37795	71.0 = 2.79527
.8 = .03149	do.0 = 1.41732	72.0 = 2.83464
.9 = .03543	37.0 = 1.45669	73.0 = 2.87401
1.0 = .03937	38.0 = 1.49606	74.0 = 2.91De8
2.0 = .07874	38.1 = 1.5 in	75.0 = 2.95275
3.0 = .11811	39.0 = 1.53543	76.0 = 2.99212
4.0 = .15748	40.0 = 1.57480	76.2 = 3.0 in
5.0 = .19685	41.0 = 1.61417	77.0 = 3.03149
6.0 = .2do22	42.0 = 1.65354	78.0 = 3.07056
7.0 = .27559	43.0 = 1.69291	79.0 = 3.11023
8.0 = .31496	44.0 = 1.73228	80.0 = 3.14960
9.0 = .354De	45.0 = 1.77165	81.0 = 3.18897
10.0 = .39370	46.0 = 1.81102	82.0 = 3.22834
11.0 = .4De07	47.0 = 1.85039	83.0 = 3.26771
12.0 = .47244	48.0 = 1.88976	84.0 = 3.30708
12.7 = .5 in	49.0 = 1.92913	85.0 = 3.34645
13.0 = .51181	50.0 = 1.96350	86.0 = 3.38582
14.0 = .55118	50.8 = 2.0 in	87.0 = 3.42519
15.0 = .59055	51.0 = 2.00787	88.0 = 3.46456
16.0 = .62992	52.0 = 2.04724	88.9 = 3.5 in
17.0 = .66929	53.0 = 2.08661	89.0 = 3.50393
18.0 = .70866	54.0 = 2.12598	90.0 = 3.54De0
19.0 = .74803	55.0 = 2.16535	91.0 = 3.58267
20.0 = .78740	56.0 = 2.20472	92.0 = 3.62204
21.0 = .82677	57.0 = 2.24409	93.0 = 3.66141
22.0 = .86614	58.0 = 2.28346	94.0 = 3.70078
23.0 = .90551	59.0 = 2.32283	95.0 = 3.74015
24.0 = .94488	60.0 = 2.do220	96.0 = 3.77952
25.0 = .98425	61.0 = 2.40157	97.0 = 3.81889
25.4 = 1.0 in	62.0 = 2.44094	98.0 = 3.85826
26.0 = 1.02do2	63.0 = 2.48031	99.0 = 3.89763
27.0 = 1.06299	63.5 = 2.5 in	100.0 = 3.93700
28.0 = 1.102do	64.0 = 2.51968	101.6 = 4.0 in







Needle Tubing Equivalents

SIZE GAUGE NO.	APPROX. O.D.	
	MM	INCHES
26	0.4	.018
24	0.5	.022
22	0.7	.028
20	0.9	.035
18	1.3	.053
16	1.6	.065
14	2.1	.083
12	2.8	.109
10	3.4	.134
8	4.2	.165
6	5.1	.203

Suction Tubing Equivalents

SIZE FRENCH GAUGE NO.	APPROX. O.D.	
	MM	INCHES
3	1.0	.039
4	1.3	.053
5	1.7	.066
6	2.0	.079
7	2.3	.092
8	2.7	.105
9	3.0	.118
10	3.3	.131
11	3.7	.144
12	4.0	.158
13	4.3	.170
14	4.7	.184
15	5.0	.197
16	5.3	.210

TABLE OF CONTENTS

De	Dental Diagnostic Instruments Zahnärztliche Diagnostikinstrumente Instrumentos De Dentales Diagnóstico Instruments De Dentaire Diagnostic Dentale Diagnostica Strumenti		De	Endodontic Instruments Endodontie-instrumente Instrumentos Endodónticos Instruments D'endodontie Strumenti Endodontici	
• Page • Seite • página • page • pagina 1 - 18		• Page • Seite • página • page • pagina 107 - 123			
37	Dental Extraction Instruments Zahnärztliche Instruments Dentales Extracción Instrumentos Instruments De Dentaire Extraction Dentale Estrazione Strumenti		38	Dental Restorative Instruments Zahnärztliche Aurative Instrumente Dentales Restaurativa Instrumentos Instruments De Dentaire Restauration Dentale Restorative Strumenti	
• Page • Seite • página • page • pagina 19 - 83		• Page • Seite • página • page • pagina 125 - 185			
37	Periodontal Instruments Parodontitis Instrumente Periodontales Instrumentos Instruments De Parodontales Strumenti Parodontali		do	Orthodontic Instruments Kieferorthopädische Instrumente Instrumentos De Ortodoncia Orthodontie Instruments Strumenti Ortodontici	
• Page • Seite • página • page • pagina 85 - 105		• Page • Seite • página • page • pagina 187 - 200			



DENTAL DIAGNOSTIC INSTRUMENTS

Sunmark Dental Diagnostic Instruments have a proven track record of quality and durability, in even the most Demanding practice. Our range covers explorers, probes, mirrors and dressing pliers with a complete selection of industry standard patterns and some unique ergonomic Designs that will improve the diagnosis procedure, provide accurate out comes and Deliver long-lasting return on investment.

Zahnärztliche Diagnostikinstrumente

Instrumentos De Dentales Diagnóstico

Instruments De Dentaire Diagnostic

Dentale Diagnostica Strumenti

Explorers - Sickle Probes
 ZahnsonDen
 Sondas exploradoras
 SonDes exploratrices
 Esploratori - SonDe falce



Explorer Single EnDed
 Fig. 3

Catalog #	Size
33-011-03	Standard



Explorer Single EnDed
 Fig. 3A

Catalog #	Size
33-013-03	Standard



Explorer Single EnDed
 Fig. 6

Catalog #	Size
33-015-01	Standard



Explorer Single EnDed
 Fig. 6A

Catalog #	Size
33-015-02	Standard



Explorer Single EnDed
 Fig. 6L

Catalog #	Size
33-015-03	Standard



Explorers - Sickle Probes
 ZahnsonDen
 Sondas exploradoras
 SonDes exploratrices
 Esploratori - SonDe falce

Explorer Single EnDed

Fig. 6XL

Catalog #	Size
33-015-04	Standard


Explorer Single EnDed

Fig. 8

Catalog #	Size
33-017-08	Standard


Explorer Single EnDed

Fig. 9

Catalog #	Size
33-019-09	Standard



Shepherd's Hook Explorer Single EnDed

Fig. 17

Catalog #	Size
33-021-17	Standard



Explorer Single EnDed

Fig. 23

Catalog #	Size
33-023-23	Standard



Explorers - Sickle Probes
 ZahnsonDen
 Sondas exploradoras
 SonDes exploratrices
 Esploratori - SonDe falce




Explorer Single EnDed
 Fig. 54

Catalog #	Size
33-025-54	Standard



Explorer double EnDed
 Fig. 1

Catalog #	Size
33-027-01	Standard



Explorer double EnDed
 Fig. 2

Catalog #	Size
33-029-02	Standard



Composite™ Explorer double EnDed
 Fig. 2



Catalog #	Size
33-029-02	Standard



Explorer double EnDed
 Fig. 2A

Catalog #	Size
33-029-13	Standard



Explorers - Sickle Probes
 ZahnsonDen
 Sondas exploradoras
 SonDes exploratrices
 Esploratori - SonDe falce

Explorer double EnDed

Fig. 2H

Catalog #	Size
33-029-26	Standard


Explorer double EnDed

Fig. 2R/2L

Catalog #	Size
33-031-22	Standard


Explorer double EnDed

Fig. 3

Catalog #	Size
33-033-03	Standard


Explorer double EnDed

Fig. 3ES

Catalog #	Size
33-035-03	Standard


Cowhorn Explorer double EnDed

Fig. 3CH

Catalog #	Size
33-037-03	Standard



Explorers - Sickle Probes
 ZahnsonDen
 Sondas exploradoras
 SonDes exploratrices
 Esploratori - SonDe falce



Composite™ Cowhorn Explorer double EnDed
 Fig. 3CH



Catalog #	Size
33-037-03 C	Standard



Explorer double EnDed
 Fig. 5


Catalog #	Size
33-039-05	Standard



Composite™ Explorer double EnDed
 Fig. 5




Catalog #	Size
33-039-05 C	Standard



Wilkins/Tufts Explorer double EnDed
 Fig. TU17/23


Catalog #	Size
33-041-17	Standard



Composite™ Wilkins/Tufts Explorer double EnDed
 Fig. TU17/23



Catalog #	Size
33-041-17 C	Standard



Explorers - Sickle Probes
 ZahnsonDen
 Sondas exploradoras
 SonDes exploratrices
 Esploratori - SonDe falce

Explorer double EnDed

Fig. 6/3

Catalog #	Size
33-043-63	Standard



Explorer double EnDed

Fig. 23/6

Catalog #	Size
33-045-23	Standard



Explorer double EnDed

Fig. EN1

Catalog #	Size
33-047-01	Standard



Explorer double EnDed

Fig. 3A/TU17

Catalog #	Size
33-049-31	Standard



Old dominion University Explorer double EnDed

Fig. 11/12

Catalog #	Size
33-051-11	Standard



PerioDONTAL Probes
 Parodontometer
 Sondas paroDONTALes
 SonDes paroDONTALes
 SonDe paroDONTALi



Measuring Probe Single EnDed
 Fig. CP2

Catalog #	Size
33-061-02	2-4-6-8-10-12



Measuring Probe Single EnDed
 Fig. CP-8

Catalog #	Size
33-061-08	3-6-8-11



Measuring Probe Single EnDed
 Fig. CP-10

Catalog #	Size
33-061-10	1-2-3-5-7-8-9-10



Measuring Probe Single EnDed
 Fig. CP-11

Catalog #	Size
33-061-11	3-6-8-11



Measuring Probe Single EnDed
 Fig. CP-12

Catalog #	Size
33-061-12	3-6-9-12



PerioDONTAL Probes
 Parodontometer
 Sondas paroDONTALes
 SonDes paroDONTALes
 SonDe paroDONTALi

Measuring Probe Single EnDed

Fig. CP-18



Catalog #	Size
33-061-18	3-5-8-10

University of North Carolina Measuring Probe Single EnDed

Fig. CP-12



Catalog #	Size
33-063-12	1-2-3-4-5-6-7-8-9-10-11-12

University of North Carolina Measuring Probe Single EnDed

Fig. CP-15



Catalog #	Size
33-063-15	1-2-3-4-5-6-7-8-9-10-11-12-13-14-15

Michigan Measuring Probe Single EnDed

Fig. PQ-0



Catalog #	Size
33-065-01	3-6-8

Williams Measuring Probe Single EnDed

Fig. PQ-W



Catalog #	Size
33-067-01	1-2-3-5-7-8-9-10

PerioDONTAL Probes
 Parodontometer
 Sondas paroDONTALes
 SonDes paroDONTALes
 SonDe paroDONTALi



WHO (CPITN) Screening Probe Single Ended
 Fig. CP-11.5B

Catalog #	Size
33-069-11	3.5-5.5-8.5-11.5



WHO (CPITN) Screening Probe double Ended
 Fig. UNC15

Catalog #	Size
33-071-15	1-2-3-4-5-6-7-8-9-10-11-12-13-14-15



WHO (CPITN) Screening Probe double Ended
 Fig. UNC15

Catalog #	Size
33-073-12	3-6-9-12



WHO Thin Williams (CPITN) Screening Probe double Ended

Catalog #	Size
33-075-26	1-2-3-5-7-8-9-10



Gickman Measuring Probe Single Ended
 Fig. 26 G

Catalog #	Size
33-077-01	1-2-3-5-7-8-9-10



PerioDONTAL Probes
 Parodontometer
 Sondas paroDONTALes
 SonDes paroDONTALes
 SonDe paroDONTALi

Goldman-Fox Measuring Probe Single EnDed


Catalog #	Size
33-079-01	1-2-3-5-7-8-9-10 FLAT

Merritt Measuring Probe Single EnDed

Fig. B



Catalog #	Size
33-081-01	1-2-3-5-7-8-9-10

Michigan Measuring Probe Single EnDed

Fig. "O"



Catalog #	Size
33-083-01	3-6-8

Thin Williams Measuring Probe Single EnDed


Catalog #	Size
33-085-01	1-2-3-5-7-8-9-10

Williams Measuring Probe Single EnDed


Catalog #	Size
33-087-01	1-2-3-5-7-8-9-10

PerioDONTAL Probes
 Parodontometer
 Sondas paroDONTALes
 SonDes paroDONTALes
 SonDe paroDONTALi



Goldman-Fox Williams Measuring Probe double Ended

Catalog #	Size
33-089-01	1-2-3-5-7-8-9-10 FLAT



Williams Explorer/Measuring Probe
 Fig. 23/QW

Catalog #	Size
33-091-23	10-9-8-7-5-3-2-1



Thin Williams Explorer/Measuring Probe
 Fig. 23/QOW

Catalog #	Size
33-093-23	10-9-8-7-5-3-2-1



Thin Williams Explorer/Measuring Probe
 Fig. 23/CP-8

Catalog #	Size
33-095-23	11-8-6-3



Explorer/Measuring Probe
 Fig. 23/CP-11

Catalog #	Size
33-095-11	11-8-6-3



PerioDONTAL Probes
 Parodontometer
 Sondas paroDONTALes
 SonDes paroDONTALes
 SonDe paroDONTALi

Explorer/Screening Probe
 Fig. 23/CP-11.58

Catalog #	Size
33-097-23	11.5-8.5-5.5-3.5



Explorer/Measuring Probe
 Fig. 23/CP-12

Catalog #	Size
33-097-12	12-9-6-3



University of North Carolina Explorer/Measuring Probe
 Fig. 23/UNC12

Catalog #	Size
33-099-12	12-11-10-9-8-7-6-5-4-3-2-1



University of North Carolina Explorer/Measuring Probe
 Fig. 23/UNC15

Catalog #	Size
33-099-15	15-14-13-12-11-10-9-8-7-6-5-4-3-2-1



Explorer/Measuring Probe
 Fig. TU17/UNC15

Catalog #	Size
33-101-17	15-14-13-12-11-10-9-8-7-6-5-4-3-2-1



PerioDONTAL Probes
 Parodontometer
 Sondas paroDONTALes
 SonDes paroDONTALes
 SonDe paroDONTALi



Williams Explorer/Measuring Probe
Fig. 23/W

Catalog #	Size
33-103-23	10-9-8-7-5-3-2-1



Thin Williams Explorer/Measuring Probe
Fig. 6/OW

Catalog #	Size
33-105-06	10-9-8-7-5-3-2-1



Williams Explorer/Measuring Probe
Fig. 6/W

Catalog #	Size
33-105-13	10-9-8-7-5-3-2-1



Thin Williams Explorer/Measuring Probe
Fig. 17/OW

Catalog #	Size
33-107-01	10-9-8-7-5-3-2-1



Explorer/Measuring Probe
Fig. 3A/UNC15

Catalog #	Size
33-109-15	15-14-13-12-11-10-9-8-7-6-5-4-3-2-1



PerioDONTAL Probes
 Parodontometer
 Sondas paroDONTALes
 SonDes paroDONTALes
 SonDe paroDONTALi

Goldman-Fox Explorer/Measuring Probe

Fig. 23/GF

Catalog #	Size
33-111-23	10-9-8-7-5-3-2-1 FLAT


Michigan Explorer/Measuring Probe

Fig. 23/O

Catalog #	Size
33-113-23	8-6-3


Thin Williams Explorer/Measuring Probe

Fig. 23/OW

Catalog #	Size
33-115-23	10-9-8-7-5-3-2-1


Thin Williams Explorer/Measuring Probe

Fig. 20F/QOW

Catalog #	Size
33-117-02	10-9-8-7-5-3-2-1


Perry Dressing Plier

Fig. 1

Catalog #	Size
33-121-01	Standard



Cotton and Dressing Pliers
 Zahnpinzetten
 Pinzas para vendajes
 Prêcelles pour cotton
 Pinzette per medicazione



London College Dressing Plier
 Fig. 2

Catalog #	Size
33-121-02	Standard



Dressing Plier
 Fig. 3

Catalog #	Size
33-121-03	Standard



Dressing Plier
 Fig. 17

Catalog #	Size
33-121-17	Standard



Dressing Plier
 Fig. 18

Catalog #	Size
33-121-18	Standard



Dressing Plier
 Fig. 18L (With Lock)

Catalog #	Size
33-123-18	Standard



Cotton and Dressing Pliers
 Zahnpinzetten
 Pinzas para vendajes
 Précelles pour cotton
 Pinzette per medicazione

Dressing Plier

Fig. 19



Catalog #	Size
33-125-19	Standard

Utility Pick-Up Dressing Plier


Catalog #	Size
33-127-00	Standard

Dressing Plier

Fig. 317 (Smooth Beaks)



Catalog #	Size
33-129-01	Standard

Meriam Dressing Plier

Fig. 318 (Serrated Beaks)



Catalog #	Size
33-129-02	Standard

Mouth Mirror (Rhodium Coated)

Fig. 4



Catalog #	Diameter
33-135-04	Ø 22 mm

Mirrors and handles
 Mun37piegel und Griffe
 Espejos y mangos
 Mirrors et manches
 Specchiotti e manici

Mouth Mirror (Rhodium Coated)
 Fig. 5

Catalog #	Diameter
33-135-05	Ø 24 mm



Cone Socket Mirror Handle

Catalog #	Size
33-136-01	Standard



Cone Socket Mirror Handle

Catalog #	Size
33-136-02	Standard





DENTAL EXTRACTION INSTRUMENTS

Zahnärztliche Instruments

Dentales Extracción Instrumentos

Instruments De Dentaire Extraction

Dentale Estrazione Strumenti

Dental Syringes, Root Elevators
 Dental-Spritzen, Wurzelheber
 Jeringas Dentales, Elevadores para raíces
 Seringues Dentaires, Elevateurs à raciness
 Siringhe Dentali, Leve per radici



Dental Aspirating Syringe

Catalog #	Size
37-150-01	Standard



Cook Waite Dental Aspirating Syringe

Catalog #	Size
37-152-01	Standard



Dental Aspirating Syringe

Catalog #	Size
37-154-01	Standard



Luxator Root Elevator Straight

Catalog #	Tip Size
37-160-02	2 mm



Luxator Root Elevator Straight

Catalog #	Tip Size
37-160-03	3 mm



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici

Luxator Root Elevator
 Straight

Catalog #	Tip Size
37-160-04	4 mm



Luxator Root Elevator
 Straight

Catalog #	Tip Size
37-160-05	5 mm



Luxator Root Elevator
 Curved

Catalog #	Tip Size
37-161-02	2 mm



Luxator Root Elevator
 Curved

Catalog #	Tip Size
37-161-03	3 mm



Luxator Root Elevator
 Curved

Catalog #	Tip Size
37-161-04	4 mm



Root Elevators
 Wurzelheber
 Elevadores para raices
 Elevateurs à racines
 Leve per radici



Luxator Root Elevator
 Curved

Catalog #	Size, Tip Size
37-161-05	Standard , 5 mm



Seldin Root Elevator
 Fig. 34 S

Catalog #	Size
37-162-34	Standard



Apexo Root Elevator
 Fig. 301

Catalog #	Size
37-164-30	Standard



Cryer Root Elevator
 Fig. 21-Mini

Catalog #	Size
37-165-21	Standard



Cryer Root Elevator
 Fig. 22-Mini

Catalog #	Size
37-165-22	Standard



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici

Cryer Root Elevator
 Fig. 45-Small

Catalog #	Size
37-165-45	Standard



Cryer Root Elevator
 Fig. 46-Small

Catalog #	Size
37-165-46	Standard



Cryer Root Elevator
 Fig. 25-Medium

Catalog #	Size
37-167-25	Standard



Cryer Root Elevator
 Fig. 26-Medium

Catalog #	Size
37-167-26	Standard



Cryer Root Elevator
 Fig. 31-Large

Catalog #	Size
37-167-31	Standard



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici



Cryer Root Elevator
 Fig. 32-Large

Catalog #	Size
37-167-32	Standard



Heidbrink Root Elevator
 Fig. 1

Catalog #	Size
37-168-01	Standard



Heidbrink Root Elevator
 Fig. 2

Catalog #	Size
37-168-02	Standard



Heidbrink Root Elevator
 Fig. 40

Catalog #	Size
37-168-40	Standard



Heidbrink Root Elevator
 Fig. 41

Catalog #	Size
37-168-41	Standard



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici

Seldin Root Elevator
 Fig. 1L



Catalog #	Size
37-169-01	Standard

Seldin Root Elevator
 Fig. 1R



Catalog #	Size
37-169-02	Standard

Seldin Root Elevator
 Fig. 2S




Catalog #	Size
37-170-20	Standard

Seldin Root Elevator
 Fig. 3S




Catalog #	Size
37-170-30	Standard

Seldin Root Elevator
 Fig. 34



Catalog #	Size
37-170-34	Standard

Root Elevators
 Wurzelheber
 Elevadores para raices
 Elevateurs à racines
 Leve per radici



Seldin Root Elevator
 Fig. 4L

Catalog #	Size
37-171-01	Standard



Seldin Root Elevator
 Fig. 4R

Catalog #	Size
37-171-02	Standard



Seldin Root Elevator
 Fig. 301W

Catalog #	Size
37-172-31	Standard



Seldin Root Elevator
 Fig. 304W

Catalog #	Size
37-172-34	Standard



Schmeckbier (Apexo) Root Elevator
 Fig. 4

Catalog #	Size
37-175-04	Standard



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici

Schmeckbier (Apexo) Root Elevator
 Fig. 5

Catalog #	Size
37-175-05	Standard



Miller (Apexo) Root Elevator
 Fig. 71

Catalog #	Size
37-177-71	Standard



Miller (Apexo) Root Elevator
 Fig. 72

Catalog #	Size
37-177-72	Standard



Miller (Apexo) Root Elevator
 Fig. 73

Catalog #	Size
37-177-73	Standard



Miller (Apexo) Root Elevator
 Fig. 74

Catalog #	Size
37-177-74	Standard



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici



Apexo Root Elevator
 Fig. 302

Catalog #	Size
37-179-32	Standard



Apexo Root Elevator
 Fig. 303

Catalog #	Size
37-179-33	Standard



Coupland Gouge
 Fig. 1

Catalog #	Size
37-182-01	Standard



Coupland Gouge
 Fig. 2

Catalog #	Size
37-182-02	Standard



Coupland Gouge
 Fig. 3

Catalog #	Size
37-182-03	Standard



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à raciness
 Leve per radici

Crane Root Pick / Elevator
 Fig. 8




Catalog #	Size
37-185-08	Standard

Stout Root Elevator
 Fig. 11



Catalog #	Size
37-186-01	Standard

Stout Root Elevator
 Fig. 11A



Catalog #	Size
37-186-02	Standard

MacMillan Gouge Root Elevator
 Fig. 11M



Catalog #	Size
37-188-11	Standard

MacMillan Gouge Root Elevator
 Fig. 12M



Catalog #	Size
37-188-12	Standard

Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici



MacMillan Gouge Root Elevator
 Fig. 13M

Catalog #	Size
37-188-13	Standard



Woodward Root Elevator
 Fig. 15

Catalog #	Size
37-189-15	Standard



Woodward Root Elevator
 Fig. 16

Catalog #	Size
37-189-16	Standard



Friedman Root Elevator
 Fig. 31F

Catalog #	Size
37-191-31	Standard



Friedman Root Elevator
 Fig. 32F

Catalog #	Size
37-191-32	Standard



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici

Root Elevator
 Fig. 38 (Serrated)



Catalog #	Size
37-193-38	Standard

Root Elevator
 Fig. 39 (Serrated)



Catalog #	Size
37-193-39	Standard

Bradley Root Elevator
 Fig. 44C



Catalog #	Size
37-194-44	Standard

Root Elevator
 Fig. 46



Catalog #	Size
37-195-01	Standard

Root Elevator
 Fig. 46R



Catalog #	Size
37-195-02	Standard

Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici



Curtis Root Elevator
 Fig. 51

Catalog #	Size
37-196-51	Standard



Curtis Root Elevator
 Fig. 52

Catalog #	Size
37-196-52	Standard



Root Elevator
 Fig. 80

Catalog #	Size
37-198-80	Standard



Root Elevator
 Fig. 81

Catalog #	Size
37-198-81	Standard



Root Elevator
 Fig. 81A

Catalog #	Size
37-200-01	Standard



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici

Root Elevator
 Fig. 81E



Catalog #	Size
37-200-02	Standard

Root Elevator
 Fig. 92 (Serrated)



Catalog #	Size
37-201-92	Standard

Root Elevator
 Fig. 93 (Serrated)



Catalog #	Size
37-201-93	Standard

Root Elevator
 Fig. 98



Catalog #	Size
37-201-98	Standard

Root Elevator
 Fig. 190



Catalog #	Size
37-203-90	Standard

Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici



Root Elevator
 Fig. 191

Catalog #	Size
37-203-91	Standard



Cogswell Root Elevator
 Fig. A

Catalog #	Size
37-204-01	Standard



Hourigan Root Elevator
 Fig. H5 (Serrated)

Catalog #	Size
37-205-01	Standard



Hourigan Root Elevator
 Fig. H6 (Serrated)

Catalog #	Size
37-205-02	Standard



Warwick-James Root Elevator
 Straight

Catalog #	Size
37-206-01	Standard



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici

Warwick-James Root Elevator
Left

Catalog #	Size
37-206-02	Standard



Warwick-James Root Elevator
Right

Catalog #	Size
37-206-03	Standard



Winter Root Elevator
Left

Catalog #	Size
37-207-01	Standard



Winter Root Elevator
Right

Catalog #	Size
37-207-02	Standard




Winter Root Elevator
Left

Catalog #	Size
37-207-03	Standard



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à raciness
 Leve per radici



Winter Root Elevator
 Right

Catalog #	Size
37-207-04	Standard



Winter Root Elevator
 Left

Catalog #	Size
37-207-05	Standard



Winter Root Elevator
 Right

Catalog #	Size
37-207-06	Standard



Winter Root Elevator
 Left

Catalog #	Size
37-207-07	Standard



Winter Root Elevator
 Right

Catalog #	Size
37-207-08	Standard



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici

Howard Root Elevator

Fig. 1



Catalog #	Size
37-209-01	Standard

Howard Root Elevator

Fig. 2



Catalog #	Size
37-209-02	Standard

Howard Root Elevator


Fig. 3



Catalog #	Size
37-209-03	Standard

Howard Root Elevator

Fig. 4



Catalog #	Size
37-209-04	Standard

Howard Root Elevator

Fig. 5



Catalog #	Size
37-209-05	Standard

Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à raciness
 Leve per radici



Howard Root Elevator
 Fig. 6

Catalog #	Size
37-209-06	Standard



Apical Root Tip Pick
 Straight

Catalog #	Size
37-210-01	Standard



Apical Root Tip Pick
 Left

Catalog #	Size
37-210-02	Standard



Apical Root Tip Pick
 Right

Catalog #	Size
37-210-03	Standard



Heidbrink Root Tip Pick
 Straight

Catalog #	Size
37-212-01	Standard



Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à racines
 Leve per radici

Heidbrink Root Tip Pick
Left



Catalog #	Size
37-212-02	Standard



Heidbrink Root Tip Pick
Right



Catalog #	Size
37-212-03	Standard




West Root Tip Pick
Straight



Catalog #	Size
37-214-01	Standard

West Root Tip Pick
Left



Catalog #	Size
37-214-02	Standard

West Root Tip Pick
Right



Catalog #	Size
37-214-03	Standard

Root Elevators
 Wurzelheber
 Elevadores para raíces
 Elevateurs à raciness
 Leve per radici



West Root Tip Pick
 Straight



Catalog #	Size
37-216-01	Standard

West Root Tip Pick
 Left




Catalog #	Size
37-216-02	Standard

West Root Tip Pick
 Right




Catalog #	Size
37-216-03	Standard

English Pattern Tooth Extraction Forcep
 Fig. 1 (For Upper Centrals, Canines and Roots)



Catalog #	Size
37-250-01	Standard

English Pattern Tooth Extraction Forcep
 Fig. 2 (For Upper Laterals)



Catalog #	Size
37-250-02	Standard

Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione

English Pattern Tooth Extraction Forcep
 Fig. 7 (For Upper Premolars and Roots)

XTRAC™

Catalog #	Size
37-250-07	Standard



English Pattern Tooth Extraction Forcep
 Fig. 13 (For Lower Premolars)

XTRAC™

Catalog #	Size
37-250-13	Standard



English Pattern Tooth Extraction Forcep
 Fig. 17 (For Upper Right Molars)

XTRAC™

Catalog #	Size
37-250-17	Standard



English Pattern Tooth Extraction Forcep
 Fig. 18 (For Upper Left Molars)

XTRAC™

Catalog #	Size
37-250-18	Standard



English Pattern Tooth Extraction Forcep
 Fig. 19 (For Upper Wisdoms)

XTRAC™


Catalog #	Size
37-250-19	Standard



Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione




Hawk's Bill English Pattern Tooth Extraction Forcep
 Fig. 22 (For Lower Molars)




Catalog #	Size
37-250-22	Standard

English Pattern Tooth Extraction Forcep
 Fig. 29 (For Upper Roots)




Catalog #	Size
37-250-29	Standard

English Pattern Tooth Extraction Forcep
 Fig. 29S (For Upper Roots, Small Beaks)




Catalog #	Size
37-251-29	Standard

English Pattern Tooth Extraction Forcep
 Fig. 30 (For Upper Roots)



Catalog #	Size
37-253-30	Standard

English Pattern Tooth Extraction Forcep (Child)
 Fig. 37 (For Upper Incisors and Canines)



Catalog #	Size
37-253-37	Standard

Extraction Forceps
Zahnzangen
Pinzas da extracción
Daviers
Pinze per estrazione

English Pattern Tooth Extraction Forcep (Child)

Fig. 39 (For Upper Molars)



Catalog #	Size
37-253-39	Standard

English Pattern Tooth Extraction Forcep

Fig. 44N (For Upper Roots)



Catalog #	Size
37-253-44	Standard

English Pattern Tooth Extraction Forcep

Fig. 51 (For Upper Wisdoms)



Catalog #	Size
37-253-51	Standard

English Pattern Tooth Extraction Forcep

Fig. 67 (For Upper Wisdoms)



Catalog #	Size
37-253-67	Standard

Hawk's Bill English Pattern Tooth Extraction Forcep

Fig. 73 (For Lower Molars)



Catalog #	Size
37-253-73	Standard

Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione



Hawk's Bill English Pattern Tooth Extraction Forcep
 Fig. 73S (For Lower Molars, Small Beaks)



Catalog #	Size
37-255-73	Standard



English Pattern Tooth Extraction Forcep
 Fig. 74 (For Lower Anteriors and Roots)



Catalog #	Size
37-255-74	Standard



English Pattern Tooth Extraction Forcep
 Fig. 74N (For Lower Anteriors and Roots, Narrow Beaks)



Catalog #	Size
37-257-74	Standard



English Pattern Tooth Extraction Forcep
 Fig. 75 (For Lower Premolars)



Catalog #	Size
37-257-75	Standard



English Pattern Tooth Extraction Forcep
 Fig. 76 (For Upper Roots)



Catalog #	Size
37-257-76	Standard



Extraction Forceps
Zahnzangen
Pinzas da extracción
Daviers
Pinze per estrazione

DENTAL EXTRACTION INSTRUMENTS

English Pattern Tooth Extraction Forcep
Fig. 76N (For Upper Roots, Narrow Beaks)



Catalog #	Size
37-259-76	Standard



English Pattern Tooth Extraction Forcep
Fig. 76S (For Upper Roots, Small Beaks)



Catalog #	Size
37-261-76	Standard



English Pattern Tooth Extraction Forcep
Fig. 79 (For Lower Wisdoms)



Catalog #	Size
37-261-79	Standard



Cowhorn English Pattern Tooth Extraction Forcep
Fig. 86 (For Lower Wisdoms)



Catalog #	Size
37-261-86	Standard



Cowhorn English Pattern Tooth Extraction Forcep
Fig. 87 (For Lower Wisdoms)



Catalog #	Size
37-261-87	Standard



Extraction Forceps
 Zahnzangen
 Pinzas de extracción
 Daviers
 Pinze per estrazione



Cowhorn English Pattern Tooth Extraction Forcep
 Fig. 89 (For Upper Right Molars)




Catalog #	Size
37-261-89	Standard



Cowhorn English Pattern Tooth Extraction Forcep
 Fig. 90 (For Upper Left Molars)



Catalog #	Size
37-261-90	Standard



English Pattern Tooth Extraction Forcep
 Fig. 94 (For Upper Right Molars)



Catalog #	Size
37-261-94	Standard



English Pattern Tooth Extraction Forcep
 Fig. 95 (For Upper Left Molars)



Catalog #	Size
37-261-95	Standard



English Pattern Tooth Extraction Forcep
 Fig. 101 (For Upper Molars and Roots)



Catalog #	Size
37-263-01	Standard



Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione

English Pattern Tooth Extraction Forcep


Fig. 107 (For Upper Canines)




Catalog #	Size
37-263-07	Standard

English Pattern Tooth Extraction Forcep


Fig. 113 (For Upper Roots)

Catalog #	Size
37-263-13	Standard

English Pattern Tooth Extraction Forcep (Child)

Fig. 123 (For Lower Anteriors, Premolars and Roots)

Catalog #	Size
37-263-23	Standard

Guy's English Pattern Tooth Extraction Forcep


Fig. 1do (For Upper Premolars and Roots)




Catalog #	Size
37-263-do	Standard

Guy's English Pattern Tooth Extraction Forcep

Fig. 137 (For Lower Centrals and Roots)

Catalog #	Size
37-263-37	Standard


Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione



English Pattern Tooth Extraction Forcep (Child)
 Fig. 138 (For Upper Anteriors, Premolars and Roots)




Catalog #	Size
37-263-38	Standard



Gibb English Pattern Tooth Extraction Forcep
 Fig. 147 (For Upper Roots)




Catalog #	Size
37-263-47	Standard



English Pattern Tooth Extraction Forcep (Child)
 Fig. 157 (For Upper Molars)




Catalog #	Size
37-263-57	Standard



English Pattern Tooth Extraction Forcep (Child)
 Fig. 158 (For Upper Molars)




Catalog #	Size
37-263-58	Standard



English Pattern Tooth Extraction Forcep (Child)
 Fig. 159 (For Upper Premolars)



Catalog #	Size
37-263-59	Standard



Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione

Hawk's Bill English Pattern Tooth Extraction Forcep (Child)

Fig. 160 (For Lower Molars)



Catalog #	Size
37-263-60	Standard



Hawk's Bill English Pattern Tooth Extraction Forcep (Child)

Fig. 161 (For Lower Molars)



Catalog #	Size
37-263-61	Standard



English Pattern Tooth Extraction Forcep (Child)

Fig. 162 (For Lower Anteriors and Roots)



Catalog #	Size
37-263-62	Standard



English Pattern Tooth Extraction Forcep (Child)

Fig. 163 (For Upper Anteriors)



Catalog #	Size
37-263-63	Standard



American Pattern Tooth Extraction Forcep

Fig. 1 (For Upper Incisors and Canines)



Catalog #	Size
37-270-01	Standard



Extraction Forceps
 Zahnzangen
 Pinzas de extracción
 Daviers
 Pinze per estrazione



Hanahan American Pattern Tooth Extraction Forcep
 Fig. 1A (For Upper Incisors and Canines)




Catalog #	Size
37-272-01	Standard



American Pattern Tooth Extraction Forcep (Child)
 Fig. 1S (For Upper Incisors and Canines)



Catalog #	Size
37-274-01	Standard



Woodward American Pattern Tooth Extraction Forcep
 Fig. 3FS (For Lower Molars)



Catalog #	Size
37-275-03	Standard



American Pattern Tooth Extraction Forcep
 Fig. 10S (For Upper Molars)



Catalog #	Size
37-275-10	Standard



American Pattern Tooth Extraction Forcep
 Fig. 16 (For Lower Molars)



Catalog #	Size
37-275-16	Standard



Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione

American Pattern Tooth Extraction Forcep (Child)

Fig. 16S (For Lower Molars)



Catalog #	Size
37-277-16	Standard


American Pattern Tooth Extraction Forcep

Fig. 17 (For Lower Molars)



Catalog #	Size
37-277-17	Standard


Harris American Pattern Tooth Extraction Forcep

Fig. 18R (For Upper Right Molars)



Catalog #	Size
37-279-17	Standard


Harris American Pattern Tooth Extraction Forcep

Fig. 18L (For Upper Left Molars)



Catalog #	Size
37-281-17	Standard


Cowhorn American Pattern Tooth Extraction Forcep

Fig. 23 (For Lower Molars)



Catalog #	Size
37-281-23	Standard



Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione



Cowhorn American Pattern Tooth Extraction Forcep (Child)
 Fig. 23S (For Lower Molars)




Catalog #	Size
37-283-23	Standard



Parmly Alveolar American Pattern Tooth Extraction Forcep
 Fig. 32 (Upper Canines, Premolars and Molars)




Catalog #	Size
37-285-32	Standard



Parmly Alveolar American Pattern Tooth Extraction Forcep
 Fig. 32A (Upper Canines, Premolars and Molars, Narrow Beaks)




Catalog #	Size
37-287-32	Standard



American Pattern Tooth Extraction Forcep (Child)
 Fig. 39 (For Upper Primary Molars)



Catalog #	Size
37-287-39	Standard



American Pattern Tooth Extraction Forcep (Child)
 Fig. 40 (For Lower Primary Molars)



Catalog #	Size
37-287-40	Standard



Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione

American Pattern Tooth Extraction Forcep
 Fig. 53R (For Upper Right Molars)



Catalog #	Size
37-287-53	Standard



American Pattern Tooth Extraction Forcep
 Fig. 53L (For Upper Left Molars)



Catalog #	Size
37-289-53	Standard



American Pattern Tooth Extraction Forcep
 Fig. 62 (For Upper and Lower Incisors and Canines)




Catalog #	Size
37-289-62	Standard



American Pattern Tooth Extraction Forcep
 Fig. 65 (For Upper Roots)



Catalog #	Size
37-289-65	Standard



Tomes American Pattern Tooth Extraction Forcep
 Fig. 69 (For Upper and Lower Roots)



Catalog #	Size
37-289-69	Standard



Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione



Nevius American Pattern Tooth Extraction Forcep
 Fig. 88R (For Upper Right Molars)



Catalog #	Size
37-291-88	Standard



Nevius American Pattern Tooth Extraction Forcep
 Fig. 88L (For Upper Left Molars)




Catalog #	Size
37-293-88	Standard



Kells American Pattern Tooth Extraction Forcep
 Fig. 99C (For Upper Incisors, Canines and Premolars)



Catalog #	Size
37-293-99	Standard



Hull American Pattern Tooth Extraction Forcep
 Fig. 101 (For Upper and Lower Premolars)



Catalog #	Size
37-295-01	Standard



Cryer American Pattern Tooth Extraction Forcep
 Fig. 150 (For Upper Incisors, Canines, Premolars and Roots)



Catalog #	Size
37-295-50	Standard



Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione

Cryer American Pattern Tooth Extraction Forcep

Fig. 150A (For Upper Incisors, Canines, Premolars, Parallel Beaks)



Catalog #	Size
37-297-50	Standard


American Pattern Tooth Extraction Forcep

Fig. 150AS (For Upper Incisors, Canines, Premolars, Split Beaks)



Catalog #	Size
37-299-50	Standard



American Pattern Tooth Extraction Forcep (Child)

Fig. 150S (For Upper Primary Teeth and Roots)



Catalog #	Size
37-301-50	Standard



Cryer American Pattern Tooth Extraction Forcep

Fig. 151 (For Lower Incisors, Canines, Premolars and Roots)



Catalog #	Size
37-301-51	Standard


Cryer American Pattern Tooth Extraction Forcep

Fig. 151A (For Lower Incisors, Canines, Premolars and Roots, Parallel Beaks)



Catalog #	Size
37-303-51	Standard



Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione




American Pattern Tooth Extraction Forcep

Fig. 151AS (For Lower Incisors, Canines, Premolars and Roots, Split Beaks)



Catalog #	Size
37-305-51	Standard



American Pattern Tooth Extraction Forcep (Child)

Fig. 151S (For Lower Primary Teeth and Roots)



Catalog #	Size
37-307-51	Standard



American Pattern Tooth Extraction Forcep

Fig. 203 (For Lower Incisors, Canines and Premolars)



Catalog #	Size
37-309-03	Standard




American Pattern Tooth Extraction Forcep

Fig. 210H (For Upper Molars)



Catalog #	Size
37-309-10	Standard




American Pattern Tooth Extraction Forcep

Fig. 210S (For Upper Molars)



Catalog #	Size
37-311-10	Standard



Extraction Forceps
 Zahnzangen
 Pinzas da extracción
 Daviers
 Pinze per estrazione

American Pattern Tooth Extraction Forcep
 Fig. 217 (For Lower Molars)




Catalog #	Size
37-311-17	Standard



American Pattern Tooth Extraction Forcep
 Fig. 222 (For Lower Molars)




Catalog #	Size
37-311-22	Standard



American Pattern Tooth Extraction Forcep
 Fig. 286 (For Upper Roots, Incisors and Premolars)



Catalog #	Size
37-311-86	Standard



Buck PerioDONTAL Knife
 Fig. 3/4

Catalog #	Size
37-321-34	Standard



Buck PerioDONTAL Knife
 Fig. 5/6

Catalog #	Size
37-321-56	Standard



Gingivectomy knives
 ParoDONTALmesser
 Bisturías para la gingivectomía
 Bistouris pour la gingivectomie
 Bisturi per gengivectomia



Crane-Kaplan PeriodONTAL Knife
 Fig. CK-3

Catalog #	Size
37-323-03	Standard



Goldman-Fox PeriodONTAL Knife
 Fig. 7

Catalog #	Size
37-325-07	Standard



Goldman-Fox PeriodONTAL Knife
 Fig. 8

Catalog #	Size
37-325-08	Standard



Goldman-Fox PeriodONTAL Knife
 Fig. 9

Catalog #	Size
37-325-09	Standard



Goldman-Fox PeriodONTAL Knife
 Fig. 11

Catalog #	Size
37-325-11	Standard



Gingivectomy knives
 ParoDONTALmesser
 Bisturíes para la gingivectomía
 Bistouris pour la gingivectomie
 Bisturi per gengivectomia

Kramer-Nevins PerioDONTAL Knife
 Fig. 7



Catalog #	Size
37-327-07	Standard

Kramer-Nevins PerioDONTAL Knife
 Fig. 11



Catalog #	Size
37-327-11	Standard

Kirkland PerioDONTAL Knife
 Fig. 15K/16K



Catalog #	Size
37-329-15	Standard

Orban PerioDONTAL Knife
 Fig. 1 /2



Catalog #	Size
37-331-12	Standard

ReDey Thin PerioDONTAL Knife
 Fig. 1



Catalog #	Size
37-333-01	Standard

Gingivectomy knives
 ParoDONTALmesser
 Bisturíes para la gingivectomía
 Bistouris pour la gingivectomie
 Bisturi per gengivectomia



ReDey Thin PerioDONTAL Knife
 Fig. 2

Catalog #	Size
37-333-02	Standard



SanDers PerioDONTAL Knife
 Fig. 1 / 2

Catalog #	Size
37-335-12	Standard



Solt PerioDONTAL Knife
 Fig. 1 / 2

Catalog #	Size
37-337-12	Standard



USC Towner PerioDONTAL Knife
 Fig. 19/20

Catalog #	Size
37-339-19	Standard



Tunnelling Knife
 Fig. 1

Catalog #	Size
37-341-01	Standard



Gingivectomy knives
 ParoDONTALmesser
 Bisturíes para la gingivectomía
 Bistouris pour la gingivectomie
 Bisturi per gengivectomia

Tunnelling Knife

Fig. 2



Catalog #	Size
37-341-02	Standard

Hoexter PerioDONTAL Curette

Fig. 1



Catalog #	Size
37-343-01	Standard

Hoexter PerioDONTAL Curette

Fig. 2



Catalog #	Size
37-343-02	Standard

Hoexter PerioDONTAL Curette


Fig. 3



Catalog #	Size
37-343-03	Standard

Kirkland PerioDONTAL Curette

Fig. 8K/9K



Catalog #	Size
37-345-89	Standard

PerioDONTAL Curettes
 ParoDONTAL Küretten
 PerioDONTALes Curesas
 ParoDONTALes Curettes
 ParoDONTALi Curette

Kramer-Nevis PerioDONTAL Curette
Fig. 3/4

Catalog #	Size
37-347-34	Standard



Kramer-Nevis PerioDONTAL Curette
Fig. 4

Catalog #	Size
37-347-40	Standard



Kramer-Mesial PerioDONTAL Curette
Fig. KRA-1

Catalog #	Size
37-349-01	Standard





Kramer-Distal PerioDONTAL Curette
Fig. KRA-2

Catalog #	Size
37-349-02	Standard



Kramer PerioDONTAL Curette
Fig. KRA-3

Catalog #	Size
37-349-03	Standard



PerioDONTAL Curettes
 ParoDONTAL Küretten
 PerioDONTALes Curetas
 ParoDONTALes Curettes
 ParoDONTALi Curette

Prichard PerioDONTAL Curette
Fig. PR1/2



Catalog #	Size
37-351-12	Standard

Sugarman PerioDONTAL Curette
Fig. 5S/6S



Catalog #	Size
37-353-01	Standard

Sugarman PerioDONTAL Curette
Fig. 9S/10S



Catalog #	Size
37-353-02	Standard

Sugarman PerioDONTAL Curette
Fig. 11S/12S



Catalog #	Size
37-353-03	Standard

Buck PerioDONTAL File
Fig. 11/12



Catalog #	Size
37-355-01	Standard

PerioDONTAL Files
 ParoDONTAL Feilen
 Limas perioDONTALes
 Limes paroDONTALes
 Lime paroDONTALi



Sugarman PerioDONTAL File
 Fig. 1S/2S

Catalog #	Size
37-357-12	Standard



Sugarman PerioDONTAL File
 Fig. 3S/4S

Catalog #	Size
37-357-34	Standard




Schluger PerioDONTAL File
 Fig. 9/10

Catalog #	Size
37-359-91	Standard



Ochsenbein PerioDONTAL Chisel
 Fig. 1

Catalog #	Size
37-360-01	Standard



Ochsenbein PerioDONTAL Chisel
 Fig. 2

Catalog #	Size
37-360-02	Standard



PerioDONTAL Chisels
 Meißel für Parodontologie
 Cinceles perioDONTALes
 Burins parodontaux
 Scalpelli paroDONTALi

Ochsenbein PerioDONTAL Chisel
 Fig. 3



Catalog #	Size
37-360-03	Standard

Ochsenbein PerioDONTAL Chisel
 Fig. 4



Catalog #	Size
37-360-04	Standard

Fedi PerioDONTAL Chisel
 Fig. 1



Catalog #	Size
37-361-01	Standard

Fedi PerioDONTAL Chisel
 Fig. 2



Catalog #	Size
37-361-02	Standard

Fedi PerioDONTAL Chisel
 Fig. 3



Catalog #	Size
37-361-03	Standard

PerioDONTAL Chisels
 Meißel für Parodontologie
 Cinceles perioDONTALes
 Burins parodontaux
 Scalpelli paroDONTALi

Fedi PerioDONTAL Chisel
Fig. 4

Catalog #	Size
37-361-04	Standard



Jovanovic PerioDONTAL Chisel
Fig. 84

Catalog #	Width Size
37-362-01	1.0 mm, 4.0 mm
37-362-02	2.0 mm, 5.0 mm
37-362-03	1.0 mm, 4 mm



RhoDes Back Action Chisel
Fig. do/37

Catalog #	Size
37-363-01	Standard



Kramer-Nevins PerioDONTAL Chisel
Fig. 1 /2

Catalog #	Size
37-363-02	Standard




Solt PerioDONTAL Chisel
Fig. 3/4

Catalog #	Size
37-363-03	Standard



PerioDONTAL Chisels
 Meißel für Parodontologie
 Cinceles perioDONTALes
 Burins parodontaux
 Scalpelli paroDONTALi

PerioDONTAL Chisel
 Fig. TG



Catalog #	Size
37-364-04	Standard

PerioDONTAL Chisel
 Fig. TG-0



Catalog #	Size
37-364-05	Standard

Remington Hoe



Catalog #	Size
37-365-06	Standard

Kirkland PerioDONTAL Chisel
 Fig. 13K/TG




Catalog #	Size
37-365-07	Standard

Kirkland PerioDONTAL Chisel
 Fig. 13K/13KL



Catalog #	Size
37-365-08	Standard

PerioDONTAL Chisels
 Meißel für Parodontologie
 Cinceles perioDONTALes
 Burins parodontaux
 Scalpelli paroDONTALi



Periotom
Fig. 1

Catalog #	Size
37-366-01	Standard



Periotom
Fig. 2

Catalog #	Size
37-366-02	Standard



Periotom
Fig. 3

Catalog #	Size
37-366-03	Standard



Periotome
Fig. 4

Catalog #	Size
37-366-04	Standard



Periotome
Fig. 5

Catalog #	Size
37-366-05	Standard



Periotomes
 Periotome
 Periotomos
 Périotomes
 Periotomi

Periotom
 Fig. 6



Catalog #	Size
37-366-06	Standard

Lucas Dental Curette
 Fig. 84



Catalog #	Size
37-367-84	Standard

Lucas Dental Curette
 Fig. 85



Catalog #	Size
37-367-85	Standard

Lucas Dental Curette
 Fig. 86



Catalog #	Size
37-367-86	Standard

Lucas Dental Curette
 Fig. 87



Catalog #	Size
37-367-87	Standard

Curettes
 Küretten
 Curetas
 Curettes
 Curettes



Lucas Dental Curette
 Fig. 88

Catalog #	Size
37-367-88	Standard




Miller Dental Curette
 Fig. 8

Catalog #	Size
37-368-08	Standard



Miller Dental Curette
 Fig. 9

Catalog #	Size
37-368-09	Standard



Miller Dental Curette
 Fig. 10

Catalog #	Size
37-368-10	Standard



Miller Dental Curette
 Fig. 11

Catalog #	Size
37-368-11	Standard



Curettes, Osteotomes
 Küretten, Osteotome
 Curetas, Osteótomos
 Curettes, Ostéotomes
 Curettes, Osteotomi

Miller Dental Curette

Fig. 12



Catalog #	Size
37-368-12	Standard

Molt Dental Curette

Fig. 2/4



Catalog #	Size
37-369-24	Standard

Bone Spreading Osteotome

Straight-Tapered Convex Tip



Catalog #	Diameter
37-370-27	Ø 2.7 mm

Bone Spreading Osteotome

Straight-Tapered Convex Tip



Catalog #	Diameter
37-370-32	Ø 3.2 mm

Bone Spreading Osteotome

Straight-Tapered Convex Tip




Catalog #	Diameter
37-370-37	Ø 3.7 mm

Osteotomes
 Osteotome
 Osteótomos
 Ostéotomes
 Osteotomi



Bone Spreading Osteotome
 Straight-Tapered Convex Tip

Catalog #	Diameter
37-370-42	Ø 4.2 mm



Bone Spreading Osteotome
 Straight-Tapered Convex Tip

Catalog #	Diameter
37-370-50	Ø 5.0 mm



Bone Spreading Osteotome
 Bayonet Shaped-Tapered Convex Tip

Catalog #	Diameter
37-371-27	Ø 2.7 mm



Bone Spreading Osteotome
 Bayonet Shaped-Tapered Convex Tip

Catalog #	Diameter
37-371-32	Ø 3.2 mm



Bone Spreading Osteotome
 Bayonet Shaped-Tapered Convex Tip

Catalog #	Diameter
37-371-37	Ø 3.7 mm



Osteotomes
 Osteotome
 Osteótomos
 Ostéotomes
 Osteotomi

Bone Spreading Osteotome
 Bayonet Shaped-Tapered Convex Tip

Catalog #	Diameter
37-371-42	Ø 4.2 mm



Bone Spreading Osteotome
 Bayonet Shaped-Tapered Convex Tip

Catalog #	Diameter
37-371-50	Ø 5.0 mm



Bone Spreading Osteotome
 Straight, Tapered Concave Tip

Catalog #	Diameter
37-372-27	Ø 2.7 mm



Bone Spreading Osteotome
 Straight, Tapered Concave Tip

Catalog #	Diameter
37-372-32	Ø 3.2 mm



Bone Spreading Osteotome
 Straight, Tapered Concave Tip

Catalog #	Diameter
37-372-37	Ø 3.7 mm



Osteotomes
 Osteotome
 Osteótomos
 Ostéotomes
 Osteotomi



Bone Spreading Osteotome
 Straight, Tapered Concave Tip

Catalog #	Diameter
37-372-42	Ø 4.2 mm



Bone Spreading Osteotome
 Straight, Tapered Concave Tip

Catalog #	Diameter
37-372-50	Ø 5.0 mm



Bone Pusher
 Straight, Slightly Concave Tip

Catalog #	Diameter
37-374-27	Ø 2.7 mm



Bone Pusher
 Straight, Slightly Concave Tip

Catalog #	Diameter
37-374-32	Ø 3.2 mm



Bone Pusher
 Straight, Slightly Concave Tip

Catalog #	Diameter
37-374-37	Ø 3.7 mm



Osteotomes
 Osteotome
 Osteótomos
 Ostéotomes
 Osteotomi

Bone Pusher

Straight, Slightly Concave Tip



Catalog #	Diameter
37-374-42	Ø 4.2 mm

Bone Pusher

Straight, Slightly Concave Tip



Catalog #	Diameter
37-374-50	Ø 5.0 mm

Bone Pusher

Bayonet Shaped



Catalog #	Size
37-375-27	Ø 2.7 mm

Bone Pusher

Bayonet Shaped



Catalog #	Size
37-375-32	Ø 3.2 mm


Bone Pusher

Bayonet Shaped



Catalog #	Size
37-375-37	Ø 3.7 mm

Osteotomes, Bone Graft Instruments
 Osteotome, Knochentransplantat Instruments
 Osteótomos, Instrumentos injerto óseo
 Ostéotomes, Instruments De greffe osseuse
 Osteotomi, Strumenti innesto osseo



Bone Pusher
 Bayonet Shaped

Catalog #	Size
37-375-42	Ø 4.2 mm



Bone Pusher
 Bayonet Shaped

Catalog #	Size
37-375-50	Ø 5.0 mm



Bone Scoop

Catalog #	Tip Size
37-380-01	5.0 mm



Palti Bone Graft Scoop

Catalog #	Tip Size
37-381-01	5 x 21 mm



Bone Graft Syringe

Catalog #	Diameter
37-383-01	Ø 5.8 mm



Bone Graft Instruments
 Knochentransplantat Instruments
 Instrumentos injerto óseo
 Instruments De greffe osseuse
 Strumenti innesto osseo

Membrane Placement Instrument



Catalog #	Size
37-385-01	Standard

LaBanca Bone Graft Plugger



Catalog #	Size, Tip Size
37-387-01	Standard , 3 x 1.5 mm

Jovanovic Bone Graft Plugger

Fig. 4



Catalog #	Size
37-389-01	Standard

Bone Packer



Catalog #	Size
37-391-01	Standard

Bone Rasp



Catalog #	Size
37-393-01	Standard

Bone Graft Instruments
 Knochentransplantat Instruments
 Instrumentos injerto óseo
 Instruments De greffe osseuse
 Strumenti innesto osseo



Grisdale Bone Condensor
 double Ended

Catalog #	Diameter
37-395-01	Ø 3 mm - 4 mm



Periosteal Raspatory / Elevator

Catalog #	Size
37-399-01	Standard



Hopkins Periosteal Raspatory / Elevator

Fig. 1

Catalog #	Size
37-399-02	Standard



Hopkins Periosteal Raspatory / Elevator

Fig. 2

Catalog #	Size
37-399-03	Standard



Mead Periosteal Raspatory / Elevator

Fig. 3

Catalog #	Size
37-399-04	Standard




Raspatories - Elevators
 Raspatorien - Elevatorien
 Legras - Elevadores
 Rugines – Elevateurs
 Periostotomi - Elevatori

Periosteal Raspatory / Elevator
 Fig. 8D



Catalog #	Size
37-399-05	Standard

Howard Periosteal Raspatory / Elevator
 Fig. 9H



Catalog #	Size
37-399-06	Standard

Allen Periosteal Raspatory / Elevator
 Fig. 9A



Catalog #	Size
37-399-07	Standard

Ohl Periosteal Raspatory / Elevator
 Fig. 10



Catalog #	Size
37-399-08	Standard

Bennet Periosteal Raspatory / Elevator
 Fig. 12S



Catalog #	Size
37-399-09	Standard

Raspatories - Elevators
 Raspatorien - Elevatorien
 Legras - Elevadores
 Rugines – Elevateurs
 Periostotomi - Elevatori



Goldman-Fox Periosteal Raspatory / Elevator
 Fig. 14



Catalog #	Size
37-399-10	Standard

Hirschfield Periosteal Raspatory / Elevator
 Fig. 20



Catalog #	Size
37-399-11	Standard

Periosteal Raspatory / Elevator
 Fig. 24G



Catalog #	Size
37-399-12	Standard

Sicilia Periosteal Raspatory / Elevator
 Fig. 24G



Catalog #	Size
37-399-13	Standard

Periosteal Raspatory / Elevator
 Fig. 149



Catalog #	Size
37-399-14	Standard

Raspatories - Elevators
 Raspatorien - Elevatorien
 Legras - Elevadores
 Rugines – Elevateurs
 Periostotomi - Elevatori

Benque Periosteal Raspatory / Elevator


Catalog #	Size
37-399-15	Standard

Woo37on Periosteal Raspatory / Elevator

Fig. 1



Catalog #	Size
37-399-16	Standard

Grisdale Periosteal Raspatory / Elevator


Catalog #	Size
37-399-17	Standard

Hourigan Periosteal Raspatory / Elevator

Fig. 2



Catalog #	Size
37-399-18	Standard

K-N Periosteal Raspatory / Elevator

Fig. 1



Catalog #	Size
37-399-19	Standard

Raspatories - Elevators
 Raspatorien - Elevatorien
 Legras - Elevadores
 Rugines – Elevateurs
 Periostotomi - Elevatori



K-N Periosteal Raspatory / Elevator
 Fig. 152

Catalog #	Size
37-399-20	Standard



Prichard Periosteal Raspatory / Elevator
 Fig. PR-3

Catalog #	Size
37-399-21	Standard



Periosteal Raspatory / Elevator
 Fig. 3/9A

Catalog #	Size
37-399-22	Standard



LaBanca Periosteal Raspatory / Elevator
 Fig. 1

Catalog #	Size
37-399-23	Standard



LaBanca Periosteal Raspatory / Elevator
 Fig. 2

Catalog #	Size
37-399-24	Standard



Raspatories - Elevators
 Raspatorien - Elevatorien
 Legras - Elevadores
 Rugines – Elevateurs
 Periostotomi - Elevatori

Buser Periosteal Raspatory / Elevator


Catalog #	Size
37-399-25	Standard

Allen Periosteal Raspatory / Elevator


Catalog #	Size
37-399-26	Standard

Periosteal Raspatory / Elevator


Catalog #	Size
37-399-27	Standard


Hoaxter Periosteal Raspatory / Elevator

Fig. 4



Catalog #	Size
37-399-28	Standard

PERIODONTAL INSTRUMENTS



Parodontitis Instrumente

PeriodONTALes Instrumentos

Instruments De ParoDONTALes

Strumenti ParoDONTALi

Raspatories - Elevators, Scalers and Curettes
 Raspatorien - Elevatorien, Scaler und Küretten
 Legras - Elevadores, Raspadores y Curetas
 Rugines - Elevateurs, Détartreurs et Curettes
 Periototomi - Elevatori, Scaler e Curette



Hygienist Sickle Scaler
 Fig. H6/H7

Catalog #	Size
37-401-67	Standard



Goldman Sickle Scaler
 Fig. 6/7

Catalog #	Size
37-403-57	Standard



Hygienist-Towner Sickle Scaler
 Fig. H5/U15

Catalog #	Size
37-405-15	Standard



Hygienist-Jacquette Sickle Scaler
 Fig. H5/33

Catalog #	Size
37-407-33	Standard



Towner-Jacquette Sickle Scaler
 Fig. U15/30

Catalog #	Size
37-409-30	Standard



Scalers and Curettes
 Scaler und Küretten
 Raspadores y Curetas
 Détartreurs et Curettes
 Scaler e Curette

Towner-Jacquette Sickle Scaler
 Fig. U15/De

Catalog #	Size
37-409-33	Standard



Jacquette Sickle Scaler
 Fig. 30/33

Catalog #	Size
37-411-30	Standard



Jacquette Sickle Scaler
 Fig. 31/32

Catalog #	Size
37-411-31	Standard



Jacquette Sickle Scaler
 Fig. 34/35

Catalog #	Size
37-411-34	Standard



Jacquette Sickle Scaler
 Fig. H3/H4

Catalog #	Size
37-413-34	Standard



Scalers and Curettes
 Scaler und K retten
 Raspadores y Curetas
 D tarteurs et Curettes
 Scaler e Curette



Posterior-Jacquette Sickle Scaler
 Fig. 14/15

Catalog #	Size
37-415-15	Standard



Sickle Scaler
 Fig. 204

Catalog #	Size
37-419-01	Standard



Sickle Scaler
 Fig. 204IUFW

Catalog #	Size
37-419-02	Standard



Sickle Scaler
 Fig. 204S

Catalog #	Size
37-419-03	Standard



Sickle Scaler
 Fig. 204SD

Catalog #	Size
37-419-04	Standard



Scalers and Curettes
 Scaler und K retten
 Raspadores y Curetas
 D tarteurs et Curettes
 Scaler e Curette

Springfield Sickle Scaler

Fig. 1



Catalog #	Size
37-421-01	Standard

Springfield Sickle Scaler

Fig. 2/3



Catalog #	Size
37-421-23	Standard

Sickle Scaler

Fig. 135



Catalog #	Size
37-423-35	Standard

Remington Sickle Scaler

Fig. 3/4



Catalog #	Size
37-425-34	Standard

Sickle Scaler

Fig. 129



Catalog #	Size
37-427-29	Standard



Scalers and Curettes
 Scaler und K retten
 Raspadores y Curetas
 D tarteurs et Curettes
 Scaler e Curette



Taylor Sickle Scaler
 Fig. T2/T3

Catalog #	Size
37-429-23	Standard



Sickle Scaler
 Fig. 137

Catalog #	Size
37-431-37	Standard



Morse Sickle Scaler
 Fig. 0/00

Catalog #	Size
37-433-00	Standard



Odont-Jacquette Sickle Scaler
 Fig. ODI/De

Catalog #	Size
37-435-De	Standard



Nebraska Sickle Scaler
 Fig. 128

Catalog #	Size
37-437-28	Standard



Scalers and Curettes
 Scaler und Küretten
 Raspadores y Curetas
 Détartreurs et Curettes
 Scaler e Curette

Sickle Scaler
 Fig. USC-128

Catalog #	Size
37-439-28	Standard



Whiteside Sickle Scaler
 Fig. WH-2

Catalog #	Size
37-441-02	Standard




Younger-Good Sickle Scaler
 Fig. YG-15

Catalog #	Size
37-443-15	Standard



Sickle Scaler
 Fig. CI-1

Catalog #	Size
37-445-01	Standard



Sickle Scaler
 Fig. CI-2/3

Catalog #	Size
37-447-23	Standard



Scalers and Curettes
 Scaler und Küretten
 Raspadores y Curetas
 Détartreurs et Curettes
 Scaler e Curette



Cumine Sickle Scaler
 Fig. 152

Catalog #	Size
37-449-52	Standard



Sickle Scaler
 Fig. B-1

Catalog #	Size
37-451-01	Standard



Sickle Scaler
 Fig. B-2/3

Catalog #	Size
37-453-23	Standard



Sickle Scaler
 Fig. B-4

Catalog #	Size
37-455-04	Standard



Sickle Scaler
 Fig. B-5/6

Catalog #	Size
37-457-56	Standard



Scalers and Curettes
 Scaler und Küretten
 Raspadores y Curetas
 Détartreurs et Curettes
 Scaler e Curette

Crane-Kaplan Sickle Scaler

Fig. CK-6

Catalog #	Size
37-459-06	Standard


Crane-Kaplan Sickle Scaler

Fig. CK-6D

Catalog #	Size
37-461-06	Standard


McCall Sickle Scaler

Fig. 11/12

Catalog #	Size
37-465-01	Standard


McCall Sickle Scaler

Fig. 11A/12A

Catalog #	Size
37-465-02	Standard



McCall Sickle Scaler

Fig. 11/12S

Catalog #	Size
37-465-03	Standard



Scalers and Curettes
 Scaler und K retten
 Raspadores y Curetas
 D tarteurs et Curettes
 Scaler e Curette



McCall Curette
 Fig. 13S/14S

Catalog #	Size
37-465-04	Standard



McCall Curette
 Fig. MTC 13/14

Catalog #	Size
37-465-05	Standard



McCall Curette
 Fig. 17S/18S

Catalog #	Size
37-465-06	Standard



McCall Curette
 Fig. 19/20

Catalog #	Size
37-465-07	Standard



Goldman-Fox Sickle Curette
 Fig. GF-1

Catalog #	Size
37-467-01	Standard



Scalers and Curettes
 Scaler und K retten
 Raspadores y Curetas
 D tarteurs et Curettes
 Scaler e Curette

Goldman-Fox Curette
 Fig. GF-2



Catalog #	Size
37-467-02	Standard

Goldman-Fox Curette
 Fig. GF-3



Catalog #	Size
37-467-03	Standard

Goldman-Fox Curette
 Fig. GF-4



Catalog #	Size
37-467-04	Standard

Buccal-Lingual Sickle Scaler
 Fig. GF-5



Catalog #	Size
37-467-05	Standard

Goldman-Fox Hoe
 Fig. GF-6



Catalog #	Size
37-467-06	Standard

Scalers and Curettes
 Scaler und Kuretten
 Raspadores y Curetas
 Détartreurs et Curettes
 Scaler e Curette



Goldman-Fox Sickle Scaler
 Fig. GF-21

Catalog #	Size
37-467-07	Standard



Goldman-Fox Sickle Scaler
 Fig. GF-21S/De

Catalog #	Size
37-467-08	Standard



Cattoni Sickle Scaler
 Fig. 103/106

Catalog #	Size
37-469-01	Standard



Cattoni Hoe/Chisel
 Fig. 104/105

Catalog #	Size
37-469-02	Standard



Cattoni Sickle Scaler
 Fig. 107/108

Catalog #	Size
37-469-03	Standard



Scalers and Curettes
 Scaler und K retten
 Raspadores y Curetas
 D tarteurs et Curettes
 Scaler e Curette

Cattoni Sickle Scaler
 Fig. 107/108S



Catalog #	Size
37-469-04	Standard

Cushing Chisel Scaler
 Fig. 61



Catalog #	Size
37-471-01	Standard

Cushing Chisel Scaler
 Fig. 62



Catalog #	Size
37-471-02	Standard

Darby-Perry Chisel
 Fig. 375K/6K



Catalog #	Size
37-473-56	Standard

Darby-Perry Chisel
 Fig. 379K/11K



Catalog #	Size
37-473-91	Standard

Scalers and Curettes
 Scaler und Kùretten
 Raspadores y Curetas
 Détartreurs et Curettes
 Scaler e Curette



Langer Curette
 Fig. 1/2

Catalog #	Size
37-475-01	Standard



Langer Curette
 Fig. 3/4

Catalog #	Size
37-475-02	Standard



Langer Curette
 Fig. 5/6

Catalog #	Size
37-475-03	Standard



Langer Curette
 Fig. 17/18

Catalog #	Size
37-475-04	Standard



Indiana University Curette
 Fig. IU13/14

Catalog #	Size
37-477-34	Standard



Scalers and Curettes
 Scaler und Küretten
 Raspadores y Curetas
 Détartreurs et Curettes
 Scaler e Curette

Indiana University Curette

Fig. IU17/18

Catalog #	Size
37-477-78	Standard


Bunting Curette

Fig. 1 /2

Catalog #	Size
37-479-12	Standard


Bunting Curette

Fig. 5/6

Catalog #	Size
37-479-56	Standard


Younger-Good Curette

Fig. YG-7/8

Catalog #	Size
37-481-56	Standard


Barnhart Curette

Fig. 1 /2

Catalog #	Size
37-483-12	Standard



Scalers and Curettes
 Scaler und Kuretten
 Raspadores y Curetas
 Détartreurs et Curettes
 Scaler e Curette



Barnhart Curette
 Fig. 5/6

Catalog #	Size
37-483-56	Standard



Columbia Curette
 Fig. 2R/2L

Catalog #	Size
37-485-02	Standard



Columbia Curette
 Fig. 4R/4L

Catalog #	Size
37-485-04	Standard



Columbia Curette
 Fig. 13/14

Catalog #	Size
37-485-13	Standard



Bates Curette
 Fig. 13S/14S

Catalog #	Size
37-487-34	Standard



Scalers and Curettes
 Scaler und K reppen
 Raspadores y Curetas
 D tarteurs et Curettes
 Scaler e Curette

Rule Curette
 Fig. 3/4

Catalog #	Size
37-489-34	Standard



Hutchinson Curette
 Fig. Hu1/ 2

Catalog #	Size
37-491-12	Standard



Ratcliff Curette
 Fig. UC3/4

Catalog #	Size
37-493-34	Standard



Lomalinda Curette
 Fig. 10/11

Catalog #	Size
37-495-01	Standard




Mallery Curette
 Fig. 1 /2

Catalog #	Size
37-497-12	Standard



Scalers and Curettes
 Scaler und K retten
 Raspadores y Curetas
 D tarteurs et Curettes
 Scaler e Curette



Curette
 Fig. ME-1

Catalog #	Size
37-499-12	Standard



Curette
 Fig. 411/412

Catalog #	Size
37-501-11	Standard



Curette
 Fig. 413/414

Catalog #	Size
37-501-13	Standard



Goldman Curette
 Fig. 31/32

Catalog #	Size
37-503-31	Standard



Gracey Curette
 Fig. 1 / 2

Catalog #	Size
37-505-01	Standard



Scalers and Curettes
 Scaler und Küretten
 Raspadores y Curetas
 Détartreurs et Curettes
 Scaler e Curette

Gracey Curette

Fig. 3/4



Catalog #	Size
37-505-02	Standard

Gracey Curette

Fig. 5/6



Catalog #	Size
37-505-03	Standard

Gracey Curette

Fig. 7/8



Catalog #	Size
37-505-04	Standard

Gracey Curette

Fig. 9/10



Catalog #	Size
37-505-05	Standard

Gracey Curette

Fig. 11/12



Catalog #	Size
37-505-06	Standard

Scalers and Curettes, Pocket Markers
 Scaler und Küretten, Taschenmarkierungspinzetten
 Raspadores y Curetas, Pinzas para marcar el sulco perioDONTAL
 Déartreurs et Curettes, Pincés à marquer la poche
 Scaler e Curette, Pinzette segnatasche



Gracey Curette
 Fig. 13/14

Catalog #	Size
37-505-07	Standard



DeMarco Furcation Curette
 Fig. 1

Catalog #	Size
37-511-01	Standard



DeMarco Furcation Curette
 Fig. 2

Catalog #	Size
37-511-02	Standard



Goldman-Fox PerioDONTAL Pocket Marker
 Left, Fig. GF-1

Catalog #	Size
37-515-01	Standard



Goldman-Fox PerioDONTAL Pocket Marker
 Right, Fig. GF-1

Catalog #	Size
37-515-02	Standard




Pocket Markers
 Taschenmarkierungspinzetten
 Pinzas para marcar el sulco periodONTAL
 Pincas à marquer la poche
 Pinzette segnatasche

Orban Medicament Loop



Catalog #	Size
37-519-00	Standard





ENDODONTIC INSTRUMENTS

Endodontie-Instrumente

instrumentos endodónticos

Instruments d'endodontie

strumenti endodontici

EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endonzia



Apical Plugger/Burnisher
 Small, Right

Catalog #	Diameter
33-521-01	Ø 0.5 mm-Ø 1.4 mm



Apical Plugger/Burnisher
 Small, Left

Catalog #	Diameter
33-521-02	Ø 0.5 mm-Ø 1.4 mm



Apical Plugger/Burnisher
 Medium, Right

Catalog #	Diameter
33-521-03	Ø 0.75 mm-Ø 1.4 mm



Apical Plugger/Burnisher
 Medium, Left

Catalog #	Diameter
33-521-04	Ø 0.75 mm-Ø 1.4 mm



Apical Plugger/Burnisher
 Large, Right

Catalog #	Diameter
33-521-05	Ø 1.0 mm-Ø 1.4 mm



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endodonzia

Apical Plugger/Burnisher

Large, Left



Catalog #	Diameter
33-521-06	Ø 1.0 mm-Ø 1.4 mm

Universal Plugger/Burnisher

Small



Catalog #	Diameter
33-523-01	Ø 0.3 mm-Ø 1.4 mm

Universal Plugger/Burnisher

Medium



Catalog #	Diameter
33-523-02	Ø 0.75 mm-Ø 1.4 mm

Universal Plugger/Burnisher

Large



Catalog #	Diameter
33-523-03	Ø 0.9 mm-Ø 1.4 mm

Explorer double EnDed

Fig. DG16



Catalog #	Size
33-531-16	Standard

EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endonzia



Explorer double EnDed
 Fig. DG16/23

Catalog #	Size
33-533-16	Standard



Explorer double EnDed
 Fig. DG16/17

Catalog #	Size
33-535-16	Standard



Glick Excavator
 Fig. 2

Catalog #	Size
33-541-02	Standard



Excavator
 Fig. 31L

Catalog #	Size
33-543-01	Standard



Excavator
 Fig. 31W

Catalog #	Size
33-543-02	Standard



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endodonzia

ENDODONTIC INSTRUMENTS

Excavator
 Fig. 32L



Catalog #	Size
33-545-01	Standard

Excavator
 Fig. DeL



Catalog #	Size
33-545-02	Standard

Locking Plier
 With Horizontal and Vertical Groove



Catalog #	Size
33-551-00	Standard

ISO Sized SpreaDer
 Fig. 30



Catalog #	Diameter Tip Length
33-553-30	Ø 0.30 mm Ø 24 mm

ISO Sized SpreaDer
 Fig. 40



Catalog #	Diameter Tip Length
33-553-40	Ø 0.40 mm Ø 24 mm



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endonzia



ISO Sized SpreaDer
 Fig. 50

Catalog #	Diameter Tip Length
33-553-50	Ø 0.50 mm Ø 28 mm



ISO Sized SpreaDer
 Fig. 60

Catalog #	Diameter Tip Length
33-553-60	Ø 0.60 mm Ø 28 mm



SpreaDer
 Fig. Gp1

Catalog #	Diameter Tip Length
33-555-01	Ø 0.25 mm Ø 20 mm



Wakai SpreaDer
 Fig. 1S

Catalog #	Diameter Tip Length
33-557-01	Ø 0.30 mm Ø 22 mm



Wakai SpreaDer
 Fig. 2S

Catalog #	Diameter Tip Length
33-557-02	Ø 0.40 mm Ø 22 mm



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endodonzia

SpreaDer
 Fig. D11

Catalog #	Diameter Tip Length
33-559-11	Ø 0.40 mm Ø 23 mm



Composite™ SpreaDer
 Fig. D11



Catalog #	Diameter
33-559-11 C	Ø 0.30 mm



SpreaDer
 Fig. D11S

Catalog #	Diameter Tip Length
33-561-11	Ø 0.35 mm Ø 23 mm



Composite™ SpreaDer
 Fig. D11S



Catalog #	Diameter
33-561-11 C	Ø 0.23 mm



SpreaDer
 Fig. D11T

Catalog #	Diameter Tip Length
33-563-11	Ø 0.23 mm Ø 23 mm



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endonzia



SpreaDer
 Fig. D11TS

Catalog #	Diameter	Tip Length
33-565-11	Ø 0.20 mm	Ø 21 mm



Composite™ SpreaDer
 Fig. D11TS




Catalog #	Diameter
33-565-11 C	Ø 0.25 mm



Composite™ SpreaDer
 Fig. 25S




Catalog #	Diameter
33-567-25 C	Ø 0.25 mm



Composite™ SpreaDer
 Fig. 25



Catalog #	Diameter
33-569-25 C	Ø 0.25 mm



Composite™ SpreaDer
 Fig. 40



Catalog #	Diameter
33-571-40 C	Ø 0.40 mm



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endonzia

Composite™ SpreaDer

Fig. 40S



Catalog #	Diameter
33-573-40 C	Ø 0.40 mm



Composite™ SpreaDer

Fig. Ma57



Catalog #	Diameter
33-575-57 C	Ø 0.20 mm



Composite™ SpreaDer

Fig. 3NT



Catalog #	Diameter
33-579-03 C	Ø 0.35 mm




Heat Carrier

Fig. 3



Catalog #	Diameter, Tip Length
33-581-03	Ø 0.35 mm, 28 mm



Heat Carrier

Fig. Gp2



Catalog #	Diameter, Tip Length
33-583-02	Ø 0.25 mm, 25 mm



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endonzia



Heat Carrier
 Fig. 0

Catalog #	Diameter, Tip Length
33-585-01	Ø 0.35 mm, 25 mm



Heat Carrier
 Fig. 00

Catalog #	Diameter, Tip Length
33-585-02	Ø 0.55 mm, 25 mm



Heat Carrier
 Fig. MA57

Catalog #	Diameter, Tip Length
33-587-57	Ø 0.20 mm, 27 mm



Heat Carrier
 Fig. Gp3

Catalog #	Diameter, Tip Length
33-589-03	Ø 0.25 mm, 30 mm



ISO Sized Plugger
 Fig. 30

Catalog #	Diameter, Tip Length
33-591-30	Ø 0.30 mm, 24 mm



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endodonzia

ISO Sized Plugger

Fig. 40

Catalog #	Diameter, Tip Length
33-591-40	Ø 0.40 mm, 24 mm



ISO Sized Plugger

Fig. 50

Catalog #	Diameter, Tip Length
33-591-50	Ø 0.50 mm, 24 mm



ISO Sized Plugger

Fig. 60

Catalog #	Diameter, Tip Length
33-591-60	Ø 0.60 mm, 24 mm



Plugger double EnDed

Fig. 1/3

Catalog #	Diameter, Tip Length
33-593-01	Ø 0.40 mm, 0.45 mm, 21 mm



Plugger double EnDed

Fig. 5/7

Catalog #	Diameter, Tip Length
33-593-02	Ø 0.50 mm, 0.75 mm, 21 mm



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endonzia



Plugger double EnDed
 Fig. 9/11

Catalog #	Diameter	Tip Length
33-593-03	Ø 1.0 mm	Ø 1.5 mm Ø 21 mm



Machtou Root Canal Plugger double EnDed
 Fig. 1 /2

Catalog #	Diameter
33-595-12	Ø 0.50 mm Ø 0.60 mm



Machtou Root Canal Plugger double EnDed
 Fig. 3/4

Catalog #	Diameter
33-595-34	Ø 0.80 mm Ø 1.10 mm



Glick BlaDe Plugger double EnDed

Catalog #	Diameter
33-597-01	Ø 1.10 mm



Retro Filling Plugger double EnDed
 Fig. 1

Catalog #	Diameter
33-599-01	Ø 0.80 mm



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endodonzia

Retro Filling Plugger double EnDed

Fig. 2



Catalog #	Diameter
33-599-02	Ø 1.10 mm

Plugger Single EnDed

Fig. 8A



Catalog #	Diameter	Tip Length
33-601-01	Ø 0.40 mm	Ø 29 mm

Plugger Single EnDed

Fig. 81/2A



Catalog #	Diameter	Tip Length
33-601-02	Ø 0.45 mm	Ø 29 mm

Plugger Single EnDed

Fig. 9A



Catalog #	Diameter	Tip Length
33-601-03	Ø 0.55 mm	Ø 29 mm

Plugger Single EnDed

Fig. 91/2A



Catalog #	Diameter	Tip Length
33-601-04	Ø 0.65 mm	Ø 29 mm

EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endonzia



Plugger Single EnDed
 Fig. 10A

Catalog #	Diameter Tip Length
33-601-05	Ø 0.75 mm Ø 29 mm



Plugger Single EnDed
 Fig. 101/2A

Catalog #	Diameter Tip Length
33-601-06	Ø 0.90 mm Ø 29 mm



Plugger Single EnDed
 Fig. 11A

Catalog #	Diameter Tip Length
33-601-07	Ø 1.05 mm Ø 29 mm



Plugger Single EnDed
 Fig. 11 1/2A

Catalog #	Diameter Tip Length
33-601-08	Ø 1.15 mm Ø 29 mm



Plugger Single EnDed
 Fig. 12A

Catalog #	Diameter Tip Length
33-601-09	Ø 1.30 mm Ø 29 mm



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endodonzia

Plugger Single EnDed
 Fig. 8

Catalog #	Diameter Tip Length
33-603-01	Ø 0.40 mm Ø 22 mm



Plugger Single EnDed
 Fig. 8 1/2

Catalog #	Diameter Tip Length
33-603-02	Ø 0.45 mm Ø 22 mm



Plugger Single EnDed
 Fig. 9

Catalog #	Diameter Tip Length
33-603-03	Ø 0.55 mm Ø 22 mm



Plugger Single EnDed
 Fig. 9 1/2

Catalog #	Diameter Tip Length
33-603-04	Ø 0.65 mm Ø 22 mm



Plugger Single EnDed
 Fig. 10

Catalog #	Diameter Tip Length
33-603-05	Ø 0.75 mm Ø 22 mm



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endonzia



Plugger Single EnDed
 Fig. 101/2

Catalog #	Diameter Tip Length
33-603-06	Ø 0.90 mm Ø 22 mm



Plugger Single EnDed
 Fig. 11

Catalog #	Diameter Tip Length
33-603-07	Ø 1.05 mm Ø 22 mm



Plugger Single EnDed
 Fig. 11 1/2

Catalog #	Diameter Tip Length
33-603-08	Ø 1.15 mm Ø 22 mm



Plugger Single EnDed
 Fig. 12

Catalog #	Diameter Tip Length
33-603-09	Ø 1.30 mm Ø 22 mm



Rubber Dam Clamp Forcep
 Fig. 4

Catalog #	Size
33-610-18	18 cm (7")



EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endodonzia

Rubber Dam Clamp Forcep


Catalog #	Size
33-613-17	17 cm (6 1/2")

**Rubber Dam Frame
 Child**


Catalog #	Size
33-615-13	13 cm (5")

**Rubber Dam Frame
 Adult**


Catalog #	Size
33-615-15	15 cm (6")

Rubber Dam Punch


Catalog #	Size
33-618-16	16 cm (6 3/4")

Ainsworth Rubber Dam Punch


Catalog #	Size
33-620-16	16 cm (6 3/4")

DENTAL RESTORATIVE INSTRUMENTS



Zahnärztliche Restaurative Instrumente

Dentales Restaurativa Instrumentos

Instruments De Dentaire Restauration

Dentale Restorative Strumenti

EnDODONTIC Instruments
 Instrumente für Endodontie
 Instrumentos para Endodoncia
 Instruments pour Endodontie
 Strumenti per Endodonzia



Composite™ Contouring Instrument



Catalog #	Size
38-651-01 C	Standard



Composite™ Burnisher
 Fig. 27/29



Catalog #	Size
38-653-27 C	Standard



Composite™ Burnisher
 Small, Medium



Catalog #	Size
38-655-01 C	Standard



Composite™ Freedman (Duckhead) Burnisher



Catalog #	Size
38-657-01 C	Standard



Composite™ Freedman Burnisher



Catalog #	Size
38-659-01 C	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Composite™ Goldfogel Embrasure Sharpening Instrument


Catalog #	Size
38-661-01 C	Standard


Composite™ Goldfogel Occlusal Anatomy Instrument


Catalog #	Size
38-661-02 C	Standard


Composite™ Goldfogel Plastic Packing Instruments


Catalog #	Size
38-661-03 C	Standard


**Composite™ Plastic Placing/ConDensing Instrument
 Medium**



Catalog #	Size
38-663-02 C	Standard


**Composite™ Plastic Placing/ConDensing Instrument
 Large**


Catalog #	Size
38-663-03 C	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Composite™ Plugger
 Fig. 5A



Catalog #	Size
38-665-05 C	Standard



Composite™ Contact Forming Instrument
 Small, Medium



Catalog #	Size
38-667-01 C	Standard



Composite™ Contact Forming Instrument
 Medium, Large



Catalog #	Size
38-667-02 C	Standard



Gol37tein Plastic Filling Instrument
 Fig. 1

Catalog #	Size
38-671-01	Standard



Composite™ Gol37tein Plastic Filling Instrument
 Fig. 1



Catalog #	Size
38-671-01 C	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Gol37tein Plastic Filling Instrument
 Fig. 2



Catalog #	Size
38-671-02	Standard

Composite™ Gol37tein Plastic Filling Instrument
 Fig. 2



Catalog #	Size
38-671-02 C	Standard

Gol37tein Plastic Filling Instrument
 Fig. 3



Catalog #	Size
38-671-03	Standard

Composite™ Gol37tein Plastic Filling Instrument
 Fig. 3



Catalog #	Size
38-671-03 C	Standard

Gol37tein Plastic Filling Instrument
 Fig. 4



Catalog #	Size
38-671-04	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Composite™ Gol37tein Plastic Filling Instrument
 Fig. 4



Catalog #	Size
38-671-04 C	Standard



Gol37tein Plastic Filling Instrument
 Fig. 5

Catalog #	Size
38-671-05	Standard



Composite™ Gol37tein Plastic Filling Instrument
 Fig. 5



Catalog #	Size
38-671-05 C	Standard



Gol37tein Plastic Filling Instrument
 Fig. 6

Catalog #	Size
38-671-06	Standard



Composite™ Gol37tein Plastic Filling Instrument
 Fig. 6



Catalog #	Size
38-671-06 C	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Composite™ Gol37tein Plastic Filling Instrument
 Micro Mini



Catalog #	Size
38-673-01 C	Standard



Gol37tein Plastic Filling Instrument
 Mini 1



Catalog #	Size
38-675-01	Standard



Composite™ Gol37tein Plastic Filling Instrument
 Mini 1



Catalog #	Size
38-675-01 C	Standard



Gol37tein Plastic Filling Instrument
 Mini 3



Catalog #	Size
38-675-02	Standard



Composite™ Gol37tein Plastic Filling Instrument
 Mini 3



Catalog #	Size
38-675-02 C	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Gol37tein Plastic Filling Instrument
 Mini 4

Catalog #	Size
38-675-03	Standard



Composite™ Gol37tein Plastic Filling Instrument
 Mini 4



Catalog #	Size
38-675-03 C	Standard



Boghosian Plastic Filling Instrument
 Fig. AB1

Catalog #	Size
38-677-01	Standard



Composite™ Boghosian Plastic Filling Instrument
 Fig. AB1



Catalog #	Size
38-677-01 C	Standard



Boghosian Plastic Filling Instrument
 Fig. AB2

Catalog #	Size
38-677-02	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Composite™ Boghosian Plastic Filling Instrument

Fig. AB2



Catalog #	Size
38-677-02 C	Standard


Plastic Filling Instrument

Fig. W3

Catalog #	Size
38-679-03	Standard


Composite™ Plastic Filling Instrument

Fig. W3



Catalog #	Size
38-679-03 C	Standard


Tufts Plastic Filling Instrument

Fig. 4F

Catalog #	Size
38-681-04	Standard


Plastic Filling Instrument

Fig. A6

Catalog #	Size
38-683-06	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Composite™ Plastic Filling Instrument
 Fig. A6



Catalog #	Size
38-683-06 C	Standard



Gregg Plastic Filling Instrument
 Fig. 4/5

Catalog #	Size
38-685-45	Standard



Composite™ Gregg Plastic Filling Instrument
 Fig. 4/5



Catalog #	Size
38-685-45 C	Standard



Interproximal Carver

Catalog #	Size
38-687-01	Standard



Composite™ Interproximal Carver



Catalog #	Size
38-687-01 C	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Composite™ Interproximal Carver


Catalog #	Size
38-687-02 C	Standard


Tufts Plastic Filling Instrument

Fig. 3

Catalog #	Size
38-689-03	Standard


Plastic Filling Instrument

Fig. 8A

Catalog #	Size
38-691-08	Standard



Composite™ Plastic Filling Instrument

Fig. 8A



Catalog #	Size
38-691-08 C	Standard


Plastic Filling Instrument

Fig. 11

Catalog #	Size
38-693-11	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Plastic Filling Instrument
 Fig. 179

Catalog #	Size
38-695-79	Standard



Plastic Filling Instrument
 Fig. G

Catalog #	Size
38-697-01	Standard



Loesche Plastic Filling Instrument
 Fig. 1

Catalog #	Size
38-699-01	Standard



Trico Plastic Filling Instrument
 Fig. 1 / 2

Catalog #	Size
38-701-12	Standard



Heidmann Plastic Spatula
 Fig. 1

Catalog #	Size
38-703-01	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Darby Perry Plastic Trimmer

Fig. 37T6



Catalog #	Size
38-705-06	Standard

Plastic Trimmer

Fig. 37T6 SE



Catalog #	Size
38-707-06	Standard

Plastic Filling Instrument

Fig. W1



Catalog #	Size
38-709-01	Standard

Woo37on Plastic Filling Instrument

Fig. 1



Catalog #	Size
38-711-01	Standard


Goldfogel Cosmetic Contouring Instrument

Fig. A



Catalog #	Size
38-715-01	Standard

Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Composite™ Goldfogel Cosmetic Contouring Instrument
 Fig. A



Catalog #	Size
38-715-01 C	Standard



Goldfogel Cosmetic Contouring Instrument
 Fig. B

Catalog #	Size
38-715-02	Standard



Composite™ Goldfogel Cosmetic Contouring Instrument
 Fig. B



Catalog #	Size
38-715-02 C	Standard



Goldfogel Cosmetic Contouring Instrument
 Fig. C

Catalog #	Size
38-715-03	Standard



Composite™ Goldfogel Cosmetic Contouring Instrument
 Fig. C




Catalog #	Size
38-715-03 C	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Goldfogel Cosmetic Contouring Instrument

Fig. D



Catalog #	Size
38-715-04	Standard

Composite™ Goldfogel Cosmetic Contouring Instrument

Fig. D



Catalog #	Size
38-715-04 C	Standard

Goldfogel Cosmetic Contouring Instrument

Fig. E



Catalog #	Size
38-715-05	Standard

Composite™ Goldfogel Cosmetic Contouring Instrument

Fig. E



Catalog #	Size
38-715-05 C	Standard

Goldfogel Cosmetic Contouring Instrument

Fig. F



Catalog #	Size
38-715-06	Standard

Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Composite™ Goldfogel Cosmetic Contouring Instrument
 Fig. F



Catalog #	Size
38-715-06 C	Standard



Goldfogel Cosmetic Contouring Instrument
 Fig. G

Catalog #	Size
38-715-07	Standard



Goldfogel Cosmetic Contouring Instrument
 Fig. H

Catalog #	Size
38-715-08	Standard



Goldfogel Cosmetic Contouring Instrument
 Fig. I

Catalog #	Size
38-715-09	Standard



Boghosian Plastic Filling Instrument
 Fig. AB1

Catalog #	Size
38-717-01	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Plastic Filling Instrument


Catalog #	Size
38-719-01	Standard

Plastic Filling Instrument

Fig. 2



Catalog #	Size
38-721-02	Standard

Baldwin Plastic Filling Instrument

Fig. 49



Catalog #	Size
38-723-49	Standard

Burnisher

Fig. 21



Catalog #	Size
38-725-21	Standard

Burnisher

Fig. 21B



Catalog #	Size
38-727-21	Standard

Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Composite™ Burnisher
 Fig. 21B



Catalog #	Size
38-727-21 C	Standard



Plastic Filling Instrument
 Fig. 1

Catalog #	Size
38-729-01	Standard



Quetin Plastic Filling Instrument
 Fig. BEQ1

Catalog #	Size
38-731-01	Standard



Plastic Filling Instrument
 Fig. 1 (Long)

Catalog #	Size
38-733-01	Standard



Plastic Filling Instrument
 Fig. 3T

Catalog #	Size
38-735-03	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Woo37on Plastic Filling Instrument
 Fig. 2



Catalog #	Size
38-737-02	Standard

Woo37on Plastic Filling Instrument
 Fig. 3



Catalog #	Size
38-737-03	Standard

Loesche Plastic Filling Instrument
 Fig. 1



Catalog #	Size
38-739-01	Standard

Loesche Plastic Filling Instrument
 Fig. 2



Catalog #	Size
38-739-02	Standard

Hatchet
 Fig. 8/9H



Catalog #	Size
38-741-89	Standard

Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Hatchet
 Fig. 13/14

Catalog #	Size
38-743-01	Standard



Hatchet
 Fig. 15/16

Catalog #	Size
38-743-02	Standard



Hatchet
 Fig. 17/18

Catalog #	Size
38-743-03	Standard



Hatchet
 Fig. 51/52

Catalog #	Size
38-743-04	Standard



Hatchet
 Fig. 53/54

Catalog #	Size
38-743-05	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Hatchet
 Fig. 19



Catalog #	Size
38-745-19	Standard

Hatchet
 Fig. 14/14



Catalog #	Size
38-747-14	Standard

Hatchet
 Fig. 15/15



Catalog #	Size
38-747-15	Standard

Hatchet
 Fig. 14/14-0



Catalog #	Size
38-749-14	Standard

Hatchet
 Fig. 15/15-0



Catalog #	Size
38-749-15	Standard

Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



WeDelstae38 Chisel
 Fig. 1 /2

Catalog #	Size
38-751-12	Standard



WeDelstae38 Chisel
 Fig. 3/4

Catalog #	Size
38-751-34	Standard



WeDelstae38 Chisel
 Fig. 5/6

Catalog #	Size
38-751-56	Standard



WeDelstaedt Chisel
 Fig. 7/10

Catalog #	Size
38-753-71	Standard



WeDelstaedt Chisel
 Fig. 8/9

Catalog #	Size
38-753-89	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

WeDelstaedt Chisel

Fig. 11/12



Catalog #	Size
38-755-11	Standard

WeDelstaedt Chisel

Fig. 40/41



Catalog #	Size
38-757-40	Standard

Hoe

Fig. 20



Catalog #	Size
38-761-20	Standard

Hoe

Fig. 21



Catalog #	Size
38-761-21	Standard

Hoe

Fig. 22



Catalog #	Size
38-761-22	Standard

Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Hoe
 Fig. 23

Catalog #	Size
38-761-23	Standard



Margin Trimmer
 Fig. 26

Catalog #	Size
38-761-26	Standard



Margin Trimmer
 Fig. 27

Catalog #	Size
38-761-27	Standard



Margin Trimmer
 Fig. 28

Catalog #	Size
38-761-28	Standard



Margin Trimmer
 Fig. 29

Catalog #	Size
38-761-29	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

DENTAL RESTORATIVE INSTRUMENTS

Margin Trimmer
 Fig. 77/78



Catalog #	Size
38-763-77	Standard

Margin Trimmer
 Fig. 79/80



Catalog #	Size
38-763-79	Standard

Angle Former
 Fig. 30/31



Catalog #	Size
38-765-30	Standard

Angle Former
 Fig. 32/De



Catalog #	Size
38-765-32	Standard

Angle Former
 Fig. 34/35



Catalog #	Size
38-765-34	Standard

Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Excavator
Fig. E1

Catalog #	Size, Diameter
38-781-01	Standard, Ø 1.2 mm



Excavator
Fig. E2

Catalog #	Diameter
38-781-02	Ø 1.5 mm



Excavator
Fig. E3

Catalog #	Diameter
38-781-03	Ø 2.5 mm



Excavator
Fig. 6

Catalog #	Diameter
38-783-06	Ø 0.85 mm



Excavator
Fig. 14

Catalog #	Diameter
38-783-14	Ø 1.5 mm



Excavators
 Excavatoren
 Excavadores
 Excavateurs
 Escavatori

DENTAL RESTORATIVE INSTRUMENTS

Excavator
 Fig. 17



Catalog #	Size, Diameter
38-783-17	Standard, Ø 1.2 mm

Excavator
 Fig. 18



Catalog #	Size, Diameter
38-783-18	Standard, Ø 1.5 mm

Excavator
 Fig. 19



Catalog #	Size, Diameter
38-783-19	Standard, Ø 2.0 mm

Excavator
 Fig. 38/39



Catalog #	Size, Diameter
38-785-38	Standard, Ø 1.5 mm

Excavator
 Fig. 17W



Catalog #	Tip Size
38-787-17	1.2 mm



Excavators
Excavatoren
Excavadores
Excavateurs
Escavatori



Excavator
Fig. 18W

Catalog #	Tip Size
38-787-18	1.5 mm



Excavator
Fig. 19W

Catalog #	Tip Size
38-787-19	1.8 mm



Darby Perry Excavator
Fig. 220/221

Catalog #	Tip Size
38-789-22	1.5 mm



Excavator
Fig. do/37

Catalog #	Tip Size
38-791-do	1.5 mm



Excavator
Fig. 40/41

Catalog #	Tip Size
38-791-40	1.1 mm



Excavators
Excavatoren
Excavadores
Excavateurs
Escavatori

DENTAL RESTORATIVE INSTRUMENTS

Excavator
Fig. 44/45

Catalog #	Tip Size
38-791-44	2.0 mm



Excavator
Fig. 61/62

Catalog #	Tip Size
38-791-61	2.0 mm



Excavator
Fig. 63/64

Catalog #	Tip Size
38-791-63	1.5 mm



Excavator
Fig. 65/66

Catalog #	Tip Size
38-791-65	1.0 mm



Excavator
Fig. 125/126

Catalog #	Diameter
38-793-25	Ø 2.5 mm



Excavators
Excavatoren
Excavadores
Excavateurs
Escavatori



Excavator
Fig. 127/128

Catalog #	Diameter
38-793-27	Ø 2.0 mm



Excavator
Fig. 129/130

Catalog #	Diameter
38-793-29	Ø 1.7 mm



Excavator
Fig. 131/132

Catalog #	Diameter
38-793-31	Ø 1.4 mm



Excavator
Fig. 153/154

Catalog #	Diameter
38-793-53	Ø 1.0 mm



Gol37tein Excavator

Catalog #	Diameter
38-795-01	Ø 1.2 mm



Excavators
 Excavatoren
 Excavadores
 Excavateurs
 Escavatori

Excavator
 Fig. 243



Catalog #	Diameter
38-797-43	Ø 1.0 mm

Excavator
 Fig. 244



Catalog #	Diameter
38-797-44	Ø 1.4 mm

Excavator
 Fig. 245



Catalog #	Diameter
38-797-45	Ø 2.0 mm

Excavator
 Fig. 246



Catalog #	Diameter
38-797-46	Ø 0.8 mm

Calcium HydroxiDe Placement Instrument



Catalog #	Size
38-799-01	Standard

Excavators
 Excavatoren
 Excavadores
 Excavateurs
 Escavatori



Calcium HydroxiDe Placement Instrument



Catalog #	Size
38-799-02	Standard

Amalgam Carrier



Catalog #	Size
38-801-01	Standard

Amalgam Carrier



Catalog #	Size
38-801-02	Standard

Amalgam Carrier



Catalog #	Size
38-801-03	Standard

Ivory Amalgam Carrier Single Ended Mini



Catalog #	Diameter
38-803-01	Ø 1.5 mm

Amalgam Instruments
 Amalgam Instrumente
 Instrumentos para amalgama
 Instruments à amalgam
 Strumenti per amalgama

Ivory Amalgam Carrier Single EnDed
 Regular

Catalog #	Diameter
38-803-02	Ø 2.0 mm



Ivory Amalgam Carrier Single EnDed
 Large

Catalog #	Diameter
38-803-03	Ø 3.2 mm



Ivory Amalgam Carrier double EnDed
 Mini/Regular

Catalog #	Diameter
38-805-01	Ø 1.5 mm Ø 2.0 mm



Ivory Amalgam Carrier double EnDed
 Regular/Large

Catalog #	Diameter
38-805-02	Ø 2.0 mm Ø 2.8 mm



Ivory Amalgam Carrier double EnDed
 Regular/Jumbo

Catalog #	Diameter
38-805-03	Ø 2.0 mm Ø 3.2 mm



Amalgam Instruments
 Amalgam Instrumente
 Instrumentos para amalgama
 Instruments à amalgam
 Strumenti per amalgama



Amalgam Well

Catalog #	Size
38-810-00	Standard



Ivory Amalgam Matrix Retainer

Catalog #	Size
38-812-01	Standard



Ivory Amalgam Matrix Retainer

Catalog #	Size
38-812-02	Standard



Tofflemire Amalgam Matrix Retainer

Catalog #	Size
38-814-01	Standard



Tofflemire Amalgam Matrix Retainer Small

Catalog #	Size
38-814-02	Standard



Matrix Retainers
 Matrizenspanner
 Portamatrices
 Porte-matrices
 Portamatrici

Andrews Plugger/ConDensor

Fig. 1 /2



Catalog #	Diameter
38-821-12	Ø 1.4 mm-Ø 2.0 mm

Black Plugger/ConDensor

Fig. 1 /2



Catalog #	Diameter
38-823-12	Ø 1.5 mm-Ø 2.5 mm

Plugger/ConDensor

Fig. 1 /2 (Serrated)



Catalog #	Diameter
38-825-12	Ø 1.5 mm-Ø 2.5 mm

Hollenback Plugger/ConDensor

Fig. H1



Catalog #	Diameter
38-827-01	Ø 1.0 mm-Ø 1.65 mm

Hollenback Plugger/ConDensor

Fig. H2



Catalog #	Diameter
38-827-02	Ø 2.0 mm - Ø 2.7 mm

Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Hollenback Plugger/ConDensor
 Fig. H3

Catalog #	Diameter
38-827-03	Ø 2.5 mm - Ø 2.5 mm



Composite™ Hollenback Plugger
 Fig. H3



Catalog #	Size
38-827-03 C	Standard



Hollenback Plugger/ConDensor
 Fig. H4

Catalog #	Diameter
38-827-04	Ø 1.8 mm - Ø 2.4 mm



Marquette Plugger/ConDensor
 Fig. 0/1 (Serrated)

Catalog #	Diameter
38-829-01	Ø 1.0 mm - Ø 1.4 mm



Marquette Plugger/ConDensor
 Fig. 0/1

Catalog #	Diameter
38-829-02	Ø 1.0 mm - Ø 1.4 mm



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Mortonson Plugger/ConDensor

Fig. 2



Catalog #	Diameter
38-831-02	Ø 0.6 mm - Ø 1.2 mm

Plugger/ConDensor

Fig. 2 (Serrated)



Catalog #	Diameter
38-833-02	Ø 1.0 mm - Ø 1.5 mm

Plugger/ConDensor

Fig. SPO (Serrated)



Catalog #	Diameter
38-835-00	Ø 0.6 mm - Ø 0.76 mm

Markley Plugger/ConDensor

Fig. 1M



Catalog #	Diameter
38-837-01	Ø 0.5 mm - Ø 0.5 mm

Markley Plugger/ConDensor

Fig. 2M




Catalog #	Diameter
38-837-02	Ø 1.1 mm - Ø 0.66 mm

Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Markley Plugger/ConDensor
 Fig. 3M



Catalog #	Diameter
38-837-03	Ø 1.5 mm - Ø 2.8 mm

Oregon Plugger/ConDensor
 Fig. 1



Catalog #	Diameter
38-839-01	Ø 0.9 mm - Ø 1.2 mm

Oregon Plugger/ConDensor
 Fig. 2



Catalog #	Diameter
38-839-02	Ø 1.0 mm - Ø 1.6 mm

Oregon Plugger/ConDensor
 Fig. 3



Catalog #	Diameter
38-839-03	Ø 1.3 mm - Ø 2.2 mm

Oregon Plugger/ConDensor
 Fig. 4



Catalog #	Diameter
38-839-04	Ø 2.0 mm - Ø 3.0 mm

Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

DENTAL RESTORATIVE INSTRUMENTS

Ward Plugger/ConDensor

Fig. 1



Catalog #	Diameter
38-841-01	Ø 1.2 mm - Ø 1.5 mm

Ward Plugger/ConDensor

Fig. 2



Catalog #	Diameter
38-841-02	Ø 2.2 mm - Ø 3.0 mm

Tanner Plugger/ConDensor

Fig. 1T



Catalog #	Diameter
38-843-01	Ø 1.7 mm - Ø 2.3 mm

Composite™ Tanner Plugger

Fig. 0T



Catalog #	Size
38-843-01 C	Standard

Tanner Plugger/ConDensor

Fig. 2T



Catalog #	Diameter
38-843-02	Ø 1.0 mm Ø 1.3 mm



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Tanner Plugger/ConDensor
 Fig. 3T

Catalog #	Diameter
38-843-03	Ø 1.6 mm Ø 2.0 mm



Plugger/ConDensor
 Fig. 9/10 (Back Action)

Catalog #	Diameter
38-845-09	Ø 1.2 mm Ø 1.6 mm



Plugger/ConDensor
 Fig. 11/12 (Back Action)

Catalog #	Diameter
38-845-11	Ø 1.5 mm Ø 2.0 mm



English Pattern Plugger/ConDensor
 Fig. 153

Catalog #	Diameter
38-847-53	Ø 1.5 mm Ø 2.5 mm



English Pattern Plugger/ConDensor
 Fig. 154

Catalog #	Diameter
38-847-54	Ø 1.4 mm Ø 2.4 mm



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Carver



Catalog #	Size
38-849-01	Standard

Carver



Catalog #	Size
38-849-02	Standard

Hollenback Carver

Fig. 1/2-3



Catalog #	Size
38-851-01	Standard

Hollenback Carver

Fig. 3S



Catalog #	Size
38-851-03	Standard

Hollenback Carver

Fig. 3



Catalog #	Size
38-853-03	Standard

Amalgam Instruments
 Amalgam Instrumente
 Instrumentos para amalgama
 Instruments à amalgam
 Strumenti per amalgama



Wiland Carver
 Fig. 8

Catalog #	Size
38-855-08	Standard



Tanner Carver
 Fig. 3

Catalog #	Size
38-857-03	Standard



Tanner Carver
 Fig. 4

Catalog #	Size
38-857-04	Standard



Tanner Carver
 Fig. 5

Catalog #	Size
38-857-05	Standard



Tanner Carver
 Fig. 6

Catalog #	Size
38-857-06	Standard



Amalgam Instruments
 Amalgam Instrumente
 Instrumentos para amalgama
 Instruments à amalgam
 Strumenti per amalgama

DENTAL RESTORATIVE INSTRUMENTS

Wall Carver
 Fig. 1



Catalog #	Size
38-859-01	Standard

Wall Carver
 Fig. 3



Catalog #	Size
38-859-03	Standard

Ward Carver
 Fig. 1



Catalog #	Size
38-861-01	Standard

Ward Carver
 Fig. 2



Catalog #	Size
38-861-02	Standard

Discoid-Cleoi37 Carver
 Fig. 1 /2



Catalog #	Size
38-863-01	Standard

Amalgam Instruments
 Amalgam Instrumente
 Instrumentos para amalgama
 Instruments à amalgam
 Strumenti per amalgama



Discoïd-Cleoi37 Carver
 Fig. 4/5

Catalog #	Size
38-863-02	Standard



Discoïd-Cleoi37 Carver
 Fig. 3/6

Catalog #	Size
38-863-03	Standard



Discoïd-Cleoi37 Carver
 Fig. 89/92

Catalog #	Size
38-863-04	Standard



Discoïd-Cleoi37 Carver
 Fig. 90B

Catalog #	Size
38-863-05	Standard



Discoïd-Cleoi37 Carver
 Fig. W-A-C-D

Catalog #	Size
38-863-06	Standard



Amalgam Instruments
 Amalgam Instrumente
 Instrumentos para amalgama
 Instruments à amalgam
 Strumenti per amalgama

Andrew Carver

Fig. 1



Catalog #	Size
38-865-01	Standard

Andrew Carver

Fig. 2



Catalog #	Size
38-865-02	Standard

Andrew Carver

Fig. 3



Catalog #	Size
38-865-03	Standard

Clappison Carver

Fig. 4



Catalog #	Size
38-865-04	Standard

Frahm Carver

Fig. 2/3



Catalog #	Size
38-867-23	Standard

Amalgam Instruments
 Amalgam Instrumente
 Instrumentos para amalgama
 Instruments à amalgam
 Strumenti per amalgama



Carver
 Fig. 104

Catalog #	Size
38-869-04	Standard



Hollenback Carver
 Fig. 1

Catalog #	Size
38-871-01	Standard



Hollenback Carver
 Fig. 2

Catalog #	Size
38-871-02	Standard



Levy Carver
 Fig. 7

Catalog #	Size
38-871-07	Standard



Shoshan Carver
 Fig. SHO-A

Catalog #	Size
38-873-00	Standard



Amalgam Instruments
 Amalgam Instrumente
 Instrumentos para amalgama
 Instruments à amalgam
 Strumenti per amalgama

Sprengel Carver

Fig. 18



Catalog #	Size
38-875-18	Standard

Tharp Carver

Fig. 1



Catalog #	Size
38-877-01	Standard

Tharp Carver

Fig. 2



Catalog #	Size
38-877-02	Standard

Rhein Amalgam File

Fig. 31/32



Catalog #	Size
38-879-31	Standard

WeDelstaedt Amalgam File

Fig. 1/4



Catalog #	Size
38-881-14	Standard

Amalgam Instruments
 Amalgam Instrumente
 Instrumentos para amalgama
 Instruments à amalgam
 Strumenti per amalgama



WeDelstaedt Amalgam File
 Fig. 2/5



Catalog #	Size
38-881-25	Standard

Carving Knife
 Fig. 2S



Catalog #	Size
38-883-02	Standard

Carving Knife
 Fig. 14L



Catalog #	Size
38-885-14	Standard

Esthetic Carving Knife
 Fig. 20



Catalog #	Size
38-887-20	Standard

Esthetic Carving Knife
 Fig. 21



Catalog #	Size
38-887-21	Standard

Amalgam Instruments
 Amalgam Instrumente
 Instrumentos para amalgama
 Instruments à amalgam
 Strumenti per amalgama

Novatech Carving Knife

Fig. 12



Catalog #	Size
38-889-12	Standard

Carving Knife

Fig. 6



Catalog #	Size
38-891-06	Standard

Black Carving Knife

Fig. 7



Catalog #	Size
38-893-07	Standard

Carving Knife

Fig. 29



Catalog #	Size
38-895-29	Standard

Ball Burnisher

Fig. 18



Catalog #	Diameter
38-901-18	Ø 2.4 mm Ø 1.6 mm

Amalgam Instruments
 Amalgam Instrumente
 Instrumentos para amalgama
 Instruments à amalgam
 Strumenti per amalgama



Ball Burnisher
 Fig. 26/27S

Catalog #	Diameter
38-901-26	Ø 2.1 mm Ø 1.3 mm



Ball Burnisher
 Fig. 27/29

Catalog #	Tip Size, Diameter
38-903-27	3.1 mm, Ø 1.6 mm



Beavertail Burnisher Single EndDed
 Fig. 2

Catalog #	Tip Size
38-905-02	2.6 mm



Beavertail Burnisher double EndDed
 Fig. 2

Catalog #	Tip Size, Diameter
38-907-02	2.6 mm, Ø 1.5 mm



Beavertail Burnisher
 Fig. 2/29

Catalog #	Tip Size
38-909-29	3.1 mm, 2.6 mm



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Acorn Burnisher
 Fig. 21



Catalog #	Diameter
38-911-01	Ø 2.7 mm, Ø 1.7 mm

Acorn Burnisher
 Fig. 21B



Catalog #	Diameter
38-911-02	Ø 3.0 mm, Ø 2.7 mm

Acorn Burnisher
 Fig. PKT-3R



Catalog #	Size
38-911-03	Standard

Romerowski Burnisher



Catalog #	Size
38-913-00	Standard

Ladmore Burnisher
 Fig. 2



Catalog #	Diameter
38-915-02	Ø 1.9 mm, Ø 1.3 mm

Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione



Composite™ Ladmore Burnisher
 Fig. 2



Catalog #	Size
38-915-02 C	Standard



Ladmore Burnisher
 Fig. 3

Catalog #	Diameter
38-915-03	Ø 1.9 mm, Ø 1.3 mm



Composite™ Ladmore Burnisher
 Fig. 3



Catalog #	Size
38-915-03 C	Standard



Miller Articulating Paper Forcep
 Fig. 3

Catalog #	Size
38-917-00	Standard



Kincheloe Gingival Retractor
 Fig. Gf10

Catalog #	Size
38-919-10	Standard



Filling Instruments
 Füllungsinstrumente
 Instrumentos para obturación
 Instruments pour obturation
 Strumenti per otturazione

Goldman-Fox Gingival Retractor

Fig. 1



Catalog #	Size
38-921-01	Standard

Meinershagen Gingival Retractor

Fig. 2



Catalog #	Size
38-923-01	Standard

Meinershagen Gingival Retractor

Fig. 3



Catalog #	Size
38-923-02	Standard

Meinershagen Gingival Retractor

Fig. 4



Catalog #	Size
38-923-03	Standard

Meinershagen Gingival Retractor


Catalog #	Size
38-923-04	Standard

Restorative Instruments
 Zahnerhaltung Instrumente
 Instrumentos para reconstituyente
 Instruments à réparation
 Strumenti per ricostituente



Gingival Cord Packer
 Fig. CSI-1

Catalog #	Size
38-925-01	Standard



Gingival Cord Packer
 Fig. 113

Catalog #	Size
38-927-13	Standard



Guyer Gingival Cord Packer
 Fig. 7

Catalog #	Size
38-929-07	Standard



Gingival Cord Packer
 Fig. BN1

Catalog #	Size
38-931-01	Standard



Gingival Cord Packer
 Fig. S6

Catalog #	Size
38-933-06	Standard



Restorative Instruments
 Zahnerhaltung Instrumente
 Instrumentos para reconstituyente
 Instruments à réparatrice
 Strumenti per ricostituente

Balshi Cord Packer


Catalog #	Size
38-935-00	Standard

Yardley Cord Packer

Fig. 1



Catalog #	Size
38-937-01	Standard

Crown Remover
 Lower


Catalog #	Size
38-941-01	Standard

Crown Remover
 Upper


Catalog #	Size
38-941-02	Standard

Crown Instruments
 Kronen Instrumente
 Instrumentos para coronas
 Instruments pour couronnes
 Strumenti per corone



Crown Remover Complete
 With Hooks Ref: 38-945-01 to 38-945-04

Catalog #	Size
38-945-00	Standard



Crown Remover Hook
 Fig. 1

Catalog #	Size
38-945-01	Standard



Crown Remover Hook
 Fig. 2

Catalog #	Size
38-945-02	Standard



Crown Remover Hook
 Fig. 3

Catalog #	Size
38-945-03	Standard



Wax Spatula
 Fig. 24

Catalog #	Size
38-950-24	Standard



Instruments For Laboratory
 Instrumente für Laboratorium
 Instrumentos para laboratorio
 Instruments pour laboratoire
 Strumenti da laboratorio

Wax Spatula
Fig. 324



Catalog #	Size
38-952-24	Standard

Wax Spatula
Fig. A6



Catalog #	Size
38-954-06	Standard

Wax Spatula
Fig. 5



Catalog #	Size
38-958-05	Standard

Wax Spatula
Fig. 7



Catalog #	Size
38-960-07	Standard

Wax Spatula



Catalog #	Size
38-962-00	Standard

Instruments For Laboratory
 Instrumente für Laboratorium
 Instrumentos para laboratorio
 Instruments pour laboratoire
 Strumenti da laboratorio



Wax Spatula
 Fig. 31



Catalog #	Size
38-964-00	Standard

Laboratory Spatula
 Fig. 7



Catalog #	Size
38-968-07	Standard

Laboratory Spatula
 Fig. 8R



Catalog #	Size
38-970-08	Standard

Laboratory Spatula
 Fig. 5A



Catalog #	Size
38-972-05	Standard

LeCron Wax Carver
 Fig. 5



Catalog #	Size
38-974-05	Standard

Instruments For Laboratory
 Instrumente für Laboratorium
 Instrumentos para laboratorio
 Instruments pour laboratoire
 Strumenti da laboratorio

Roach Wax Carver


Catalog #	Size
38-976-00	Standard

Vehe Wax Carver


Catalog #	Size
38-978-00	Standard

PK Thomas Wax Instrument

Fig. 1



Catalog #	Size
38-980-01	Standard

PK Thomas Wax Instrument

Fig. 2



Catalog #	Size
38-980-02	Standard


PK Thomas Wax Instrument

Fig. 3



Catalog #	Size
38-980-03	Standard

Instruments For Laboratory
 Instrumente für Laboratorium
 Instrumentos para laboratorio
 Instruments pour laboratoire
 Strumenti da laboratorio



PK Thomas Wax Instrument
 Fig. 4

Catalog #	Size
38-980-04	Standard



PK Thomas Wax Instrument
 Fig. 5

Catalog #	Size
38-980-05	Standard



Shaw Wax Instrument
 Fig. 1

Catalog #	Size
38-982-01	Standard



Shaw Wax Instrument
 Fig. 2

Catalog #	Size
38-982-02	Standard



Shaw Carver
 Fig. 3

Catalog #	Size
38-982-03	Standard



Instruments For Laboratory
 Instrumente für Laboratorium
 Instrumentos para laboratorio
 Instruments pour laboratoire
 Strumenti da laboratorio

Shaw Burnisher
 Fig. 4




Catalog #	Size
38-982-04	Standard

Shaw Spatula
 Fig. 7



Catalog #	Size
38-982-07	Standard



ORTHODONTIC INSTRUMENTS

Kieferorthopädische Instrumente

Instrumentos De ortodoncia

Orthodontie instruments

Strumenti ortodontici

ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Orthodontia
 Strumenti da Orthodontie

Multi Use Cutter

For Soft Wire Pins and Ligature Wires upto 0.38 mm

ULTRACut™

Catalog #	Size
36-989-01	Standard


Mini Pin and Ligature Cutter

For Soft Wire Pins and Ligature Wires upto 0.30 mm

ULTRACut™

Catalog #	Size
36-989-02	Standard


Micro Mini Pin and Ligature Cutter

For Soft Wire Pins and Ligature Wires upto 0.30 mm

ULTRACut™

Catalog #	Size
36-989-03	Standard


Micro Mini Pin and Ligature Cutter

With Long Handle - For Soft Wire Pins and Ligature Wires upto 0.30 mm

ULTRACut™

Catalog #	Size
36-989-04	Standard


Micro Mini Pin and Ligature Cutter

Angled 15° - For Soft Wire Pins and Ligature Wires upto 0.30 mm

ULTRACut™

Catalog #	Size
36-990-05	Standard



ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Orthodontia
 Strumenti da Orthodontie

Universal Distal End Cutter

With Long Handle - For Cutting All Types of Wires from 0.30 mm to 0.53 mm x 0.64 mm

ULTRA.Cut™

Catalog #	Size
36-990-01	Standard


Universal Distal End Cutter

With Long Handle - For Cutting All Types of Wires from 0.30 mm to 0.53 mm x 0.64 mm

ULTRA.Cut™

Catalog #	Size
36-990-02	Standard


Flush Cut Distal End Cutter No Hold

For All Types of Wires up to a maximum of 0.53 mm x 0.64 mm

ULTRA.Cut™

Catalog #	Size
36-990-03	Standard


Wire Cutter

Angled 15° - For hard Wires up to 0.56 mm x 0.71 mm

ULTRA.Cut™

Catalog #	Size
36-990-04	Standard


Wire Cutter

Straight - For hard Wires up to 0.56 mm x 0.71 mm

ULTRA.Cut™

Catalog #	Size
36-990-05	Standard



ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Orthodontia
 Strumenti da Orthodontie



Wire Cutter

For Cutting Soft Wire Pins and Ligature Wires up to 0.30 mm

ULTRACut™

Catalog #	Size
36-990-06	Standard



Flush Cut and Hold Distal End Cutter

For Cutting All Types of Wires up to 0.51 mm

ULTRACut™

Catalog #	Size
36-990-07	Standard



Weingart Utility Plier

ULTRAGripX™

Catalog #	Size
36-990-08	Standard



**Weingart Utility Plier
Slim**

ULTRAGripX™

Catalog #	Size
36-990-09	Standard



**How Utility Plier
Straight**

Catalog #	Size
36-990-10	Standard



ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Orthodontia
 Strumenti da Orthodontie

How Utility Plier
 Curved



Catalog #	Size
36-990-11	Standard

Band Seating Plier



Catalog #	Size
36-990-12	Standard

Adhesive Removing Plier



Catalog #	Size
36-990-13	Standard

Posterior Band Removing Plier
 Long



Catalog #	Size
36-990-14	Standard

Posterior Band Removing Plier
 Short



Catalog #	Size
36-990-15	Standard

ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Orthodontia
 Strumenti da Orthodontie



Tip Back Plier

Catalog #	Size
36-990-16	Standard



Force Module Separating Plier

Catalog #	Size
36-990-17	Standard



Coon Suture Tying Plier

Catalog #	Size
36-990-18	Standard



Crown & Band Contouring Plier

Catalog #	Size
36-990-19	Standard



Bracket Removing Plier

Catalog #	Size
36-990-20	Standard



ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Orthodontia
 Strumenti da Orthodontie

Bracket Removing Plier
 Angulated

Catalog #	Size
36-990-21	Standard



Band Removing Plier

Catalog #	Size
36-990-22	Standard



Jarabak Wire Forming Plier

Catalog #	Size
36-990-23	Standard



Three Jaw Plier

Catalog #	Size
36-990-24	Standard



Ortho Finishing Step Plier

Catalog #	Size
36-990-25	Standard



ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Orthodontia
 Strumenti da Orthodontie



Omega Loop Forming Plier

Catalog #	Size
36-990-26	Standard



Bird Beak Plier

Catalog #	Size
36-990-27	Standard



Light Wire Bird Beak Plier

Catalog #	Size
36-990-28	Standard



Tweed Loop Forming Plier

Catalog #	Size
36-990-29	Standard



Arch Bending Plier

Catalog #	Size
36-990-30	Standard



ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Ortodontia
 Strumenti da Orthodontie

Tweed Plier
 Slim

Catalog #	Size
36-990-31	Standard


Rectangular Arch Bending Plier

Catalog #	Size
36-990-32	Standard


Lingual Arch Forming Plier

Catalog #	Size
36-990-De	Standard


Hollow Chop Contouring Plier

Catalog #	Size
36-990-34	Standard


O'Brien Plier

Catalog #	Size
36-990-35	Standard



ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Orthodontia
 Strumenti da Orthodontie



Mini Three Jaw Plier

Catalog #	Size
36-990-do	Standard



Marcotte Looping Plier

Catalog #	Size
36-990-37	Standard



Hollow Chop Plier

Catalog #	Size
36-990-38	Standard



Light Wire Bird Beak Plier

Catalog #	Size
36-990-39	Standard



Nance Loop Forming Plier

Catalog #	Size
36-990-40	Standard



ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Orthodontia
 Strumenti da Orthodontie

Adams Plier

Catalog #	Size
36-990-41	Standard


V-Stop Plier

Catalog #	Size
36-990-42	Standard


Hook Crimping Plier

Catalog #	Size
36-990-43	Standard


Optical Plier

Catalog #	Size
36-990-44	Standard


Loop Forming Plier

Catalog #	Size
36-990-45	Standard



ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Orthodontia
 Strumenti da Orthodontie



Bird Beak Plier With Cutter

Catalog #	Size
36-990-46	Standard



Long Tapered Bird Beak Plier

Catalog #	Size
36-990-47	Standard



Hammerhead Tie Back Plier

Catalog #	Size
36-990-48	Standard



Ligature Director

Catalog #	Size
36-993-01	Standard



Ligature Director


Catalog #	Size
36-993-02	Standard



ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Orthodontia
 Strumenti da Orthodontie

Ligature Director/Plugger


Catalog #	Size
36-993-03	Standard

Tucker Ligature Director


Catalog #	Size
36-993-04	Standard

Plugger/Ligature Picker


Catalog #	Size
36-993-05	Standard


Fillion Explorer/Ligature Director


Catalog #	Size
36-993-06	Standard

Mershon Band Pusher


Catalog #	Size
36-995-01	Standard

ORTHODONTIC Instruments
 Instrumente für Kieferorthopädie
 Instrumentos para Ortodoncia
 Instruments pour Orthodontia
 Strumenti da Orthodontie



Band Seater



Catalog #	Size
36-995-02	Standard

Band Pusher / Scaler

Fig. 1



Catalog #	Size
36-997-01	Standard

Band Pusher / Scaler

Fig. 2



Catalog #	Size
36-997-02	Standard

Direct Bond Bracket Holder



Catalog #	Size
36-998-00	Standard



SUNMARK
INSTRUMENTS

Mark Of Quality

Our Range of Best Selling Instruments

DENTAL INSTRUMENTS CATALOGUE

+92 (331) 661 9646

www.sunmarkinstruments.com

sales@sunmarkinstruments.com

Small Industrial Estate Sialkot-51310 Punjab Pakistan

